

Introduction

The primary goal of solid phase extraction (SPE) with ProElut™ is the selective extraction of the components of interest from a complex sample or much larger sample volume prior to actual analysis (e.g. HPLC, GC). As SPE works on the principle of liquid chromatography, this is achieved by using strong but reversible interactions between the analyte and surface of the stationary phase. Typical interactions are hydrophobic (the Van der Waals force), polar (hydrogen bonding, dipole-dipole interaction) or ion exchange interactions. Interaction between the stationary phase and matrix should not occur. It is thus meaningful to carry out appropriate sample pre-treatment as this emphasizes the differences in chemical properties between the substance to be analyzed and matrix components so that these are then achieved by altering the pH or the ionic strength of the sample solution.

Under these conditions, the analyte is enriched as a narrow zone on the stationary phase. Subsequent to a washing step, which serves to remove possible adsorbed sample components, the actual selective elution of the analyte takes place.



ProElut™ SPE products are packed and assembled using custom designed equipment. Every part of ProElut™ manufacturing process is carefully monitored, we only accept products that meet our high quality standards.

SPE Principles and Techniques

SPE is a chromatographic technique first developed during the mid-1980s and is increasingly used for sample pre-treatment. The main objectives of SPE are removal of interfering matrix components and selective concentration and isolation of the analytes. This is done either by retaining the substance of interest and washing off everything else or by retaining the interfering substances and eluting the product of interest. Compared to traditional liquid / liquid extraction, SPE is more rapid, uses less solvent, eliminates emulsions, and can be automated. Additionally, a sample preparation task can often be solved more specifically by using SPE, since different interactions of the analyte with the adsorbent are possible, and methods can be optimized by adjusting chromatographic conditions. SPE offers a multitude of adsorbents for polar, hydrophobic and / or ionic interactions and has been widely used in medicine, food, environmental protection, commodity inspection, cosmetics and other fields.

The most popular SPE products are: Normal Phase, Reversed Phase, Ion Exchange Phase and Mix Mode. It is important to select the most suitable product for each application and sample.

| Type | Reversed Phase | Normal Phase | Ion Exchange Phase | Mix Mode |
|-----------------|--|--|---|--|
| Separation Mode | Reserved phase separation involves a polar (usually aqueous) or moderately polar sample matrix and a non-polar stationary phase. The analyte of interest is typically mid- to non-polar. C18 is the most common reversed phase packing | Normal phase SPE typically involves a mid- to non-polar sample matrix and a polar stationary phase. The analytes are polar compounds and the bonded phases are typically NH ₂ , PSA, polar adsorbent, silica, Florisil, and alumina, etc. | Ion exchange SPE can be used for compounds that are charged when in solution (usually aqueous, but sometimes organic). The primary retention mechanism of the compound is based on the electrostatic attraction of the charged functional group on the compound to the charged group that is bonded to the silica surface. In order for a compound to be retained by ion exchange from an aqueous solution, the pH of the sample matrix must be one at which both the compound of interest and the functional group on the bonded silica are charged The most common stationary phases are SAX and SCX | A mix of ion exchange and reserved phase retention mechanisms. This can be used for hydrophobic compounds (reversed phase) and compounds that are charged in solution (ion exchange). The most common stationary phases are PXC, PXA, etc. |
| Sample | Polar or hydrophobic compounds | Polar compounds | Compounds with charges | Compounds with charges, polar or hydrophobic compounds |
| Elution | Water or organic solvents | Non-polar solvents | Water or non-polar organic solvents | Water or organic solvents |

Benefits of SPE

- Switch sample matrices to a form more compatible with chromatographic analyses
- Concentrate analytes for increased sensitivity
- Remove interference to simplify chromatography and improve quantitation
- Protect the analytical column from contaminants


Use SPE for Samples That

- Contain particulate matter that may cause system clogging and high back pressure
- Contain components that cause high background, misleading peaks, and / or poor sensitivity
- Require cleanup, trace enrichment / concentration, or purification
- Require sample matrix or solvent exchange

Choosing a SPE Product

1. *Characterize the sample.* Factors such as the analyte's polarity relative to the matrix, the presence of charged functional groups, solubility, molecular weight, etc., determine how strongly the analyte is retained by the packed bed.
2. *Select a retention strategy.* Two approaches are possible: retain interfering compounds while the analyte passes through, or retain the analyte while interfering compounds pass through. This second approach allows concentration of the sample during analyte elution.
3. *Select proper packing type and bed volume.* Choosing the packing material with the proper selectivity results in the cleanest extract with the highest recovery. Poor sample recovery often occurs when the packed bed dimensions are not optimized. Too large of a bed volume results in an incomplete elution while too small of a bed volume results in an incomplete retention. Due to the unknown composition of many samples, experimentation may be required to determine the optimum bed dimensions for an application. Start with an intermediate bed volume, such as 200 mg or 500 mg. If you observe complete retention, you may be able to use a smaller bed volume and elution volume. If you observe incomplete retention, you will need to use a larger bed volume and elution volume.
4. *Select suitable conditioning, wash, and elution solvents.* Consider the solvent strength relative to the packing material. The final conditioning solvent should be weak, so as not to act as an eluting solvent. Buffers should be used to control ionization of potentially charged compounds. Wash solvents should remove weakly retained interferences without being strong enough to elute the analyte. Elution solvents should be strong enough to completely elute an analyte in a small volume.

SPE Phase Selection



ProElut™ SPE

Common SPE Applications

- Pharmaceutical compounds and metabolites in biological fluids
- Drugs of abuse in biological fluids
- Environmental pollutants in drinking and waste water
- Pesticides and antibiotics in food / agricultural matrices
- Desalting of proteins and peptides
- Fractionation of lipids
- Water and fat soluble vitamins

General SPE Procedures

Reversed Phase

(extraction of hydrophobic or polar organic analytes from aqueous matrix)

A. Conditioning

Rinse tube with 3 - 5 mL of methanol follow by 3 - 5 mL of deionized water / buffer (do not allow tube to dry before next step).

B. Sample Application

Apply sample to the top of the tube and draw through the packing bed.

C. Tube Wash

Wash with 5 mL of a polar solvent if analyte is to be retained (deionized water, buffer or aqueous / organic mixtures are most often used).

D. Elution

Elute analyte into a collection tube with 1 - 5 mL of a non-polar solvent.

Normal Phase

(extraction of polar analytes from non-polar organic solvents)

A. Conditioning

Rinse tube with 3 - 5 mL of non-polar solvent.

B. Sample Application

Apply sample to the top of the tube and draw through the packing bed.

C. Tube Wash

Wash with 5 mL of a non-polar solvent if analyte is to be retained.

D. Elution

Elute analyte into a collection tube with 1 - 5 mL of a polar solvent.

Ion Exchange

(extraction of charged analytes from aqueous or non-polar organic samples)

A. Conditioning

Rinse tube with 3 - 5 mL of deionized water or low ionic strength buffer (10 mM).

B. Sample Application

Apply sample to the top of the tube and draw through the packing bed (ion exchange kinetics is slower than reverse or normal phase, so keep the flow slow).

C. Tube Wash

Wash with 5 mL of deionized water or low ionic strength buffer.

D. Elution

Elute analyte into a collection tube with 1 - 5 mL of buffer at high ionic strength (0.1 - 1 M) or modified pH (pH such that the analyte is uncharged).


Features of ProElut™ SPE

- Rapid sample preparation within minutes
- Higher recoveries without the formation of emulsion
- High precision of analytical results by use of disposable cartridges
- Saving of solvent and hence reduction in both materials costs and cost of disposal
- Possibilities for automating the entire process
- Optimized, validated and certified manufacturing

ProElut™ SPE Features and Benefits

ProElut™ SPE uses spherical silica with high purity and narrow particle size distribution as support. The spherical silica with fewer fines gives a more regular, stable and reproducible chromatography bed that gives a faster, more even flow rate for better separation. Fines cause back pressure increases that can result in clogging and can pass through filters and contaminate the final product. A narrower particle size distribution will give a more homogenous packing that will help in collecting more concentrated fractions and reducing solvent consumption to achieve higher recovery and reproducible results.

Figure 1. Comparison of Silica Shape and Particle Size Distribution




High Purity Silica Improves Recovery and Reproducibility

- ProElut™ packings are considerably lower in trace metals than the other materials
- High purity silica to enhance retention capacity of basic compounds
- Irregular silica, depending on its method of manufacturing, normally contains trace quantities of a variety of different metals, which in turn can affect the separation

Table 1. Metal Analysis of Base Silica (ppm)


| Metal | Dikma ProElut™ High-purity Silica | Brand A Silica | Brand B Silica |
|-------|-----------------------------------|----------------|----------------|
| Na | < 12 | 917 | 17 |
| Al | < 12 | 276 | 57 |
| Ca | < 12 | < 12 | < 12 |
| Fe | < 12 | 64 | 23 |
| Ti | < 12 | < 12 | 130 |
| Zr | < 12 | 48 | 38 |

Figure 2. Level of Organic Extractables from SPE Packing Materials as Measured by GC-FID*



*Sep-Pak is a registered trademark of Waters Corporation. Dikma Technologies Inc. is not affiliated with the above company.

Figure 3. Level of Organic Extractables of Cartridge Tubes and Frits as Measured by GC-FID*



Sample Preparation

ProElut™ SPE

ProElut™ SPE

Sorbent Specifications

| Sorbent Phase | Category | Base Material | Particle Size (µm) | Pore Size (Å) | Surface Area (m ² /g) | Bonded Functional Group | Carbon Loading | Endcapping |
|-----------------|---|-----------------|--------------------|---------------|----------------------------------|--|----------------|------------|
| PLS | Reversed phase | PS-DVB | 50 | 80 | 800 | Hydrophilic / lipophilic | - | - |
| PWC | Reversed phase / weak cation exchange | PS-DVB | 50 | 80 | 800 | Hydrophilic / lipophilic, carboxylic acid | - | - |
| PWA | Reversed phase / weak anion exchange | PS-DVB | 50 | 80 | 800 | Hydrophilic / lipophilic, ethylene diamine | - | - |
| PXC | Reversed phase / strong cation exchange | PS-DVB | 50 | 80 | 800 | Hydrophilic / lipophilic, sulfonic acid | - | - |
| PXA | Reversed phase / strong anion exchange | PS-DVB | 50 | 80 | 800 | Hydrophilic / lipophilic, quaternary amine | - | - |
| C18 | Reversed phase | Silica | 50 | 60 | 500 | Octadecyl | 19% | Yes |
| C18-U | Reversed phase | Silica | 50 | 60 | 500 | Octadecyl, silanol | 17% | No |
| C8 | Reversed phase | Silica | 50 | 60 | 500 | Octyl | 11% | Yes |
| C2 | Reversed phase | Silica | 50 | 60 | 500 | Ethyl | 5.6% | Yes |
| PH | Reversed phase | Silica | 50 | 60 | 500 | Phenyl | 10% | Yes |
| CN | Normal / reversed phase | Silica | 50 | 60 | 500 | Cyanopropyl | 8% | Yes |
| SCX | Strong cation exchange | Silica | 50 | 60 | 500 | Benzenesulfonic acid | 10.9% | No |
| SAX | Strong anion exchange | Silica | 50 | 60 | 500 | Trimethylaminopropyl | 8% | No |
| Silica | Normal phase | Silica | 50 | 60 | 500 | Silanol | - | No |
| NH ₂ | Normal phase / weak anion exchange | Silica | 50 | 60 | 500 | Aminopropyl | 5.5% | No |
| PSA | Normal phase / weak anion exchange | Silica | 50 | 60 | 500 | Ethylenediamino-N-propyl | 8.5% | No |
| AL-A | Normal phase | Acidic alumina | 125 | - | 200 | Acidic alumina | - | - |
| AL-B | Normal phase | Basic alumina | 125 | - | 200 | Basic alumina | - | - |
| AL-N | Normal phase | Neutral alumina | 125 | - | 200 | Neutral alumina | - | - |
| Florisil | Normal phase | Magnesium | 150 - 200 | - | - | - | - | - |
| CARB | Nonbonded carbon phase | Carbon | 120 - 400 | - | 100 | - | - | - |

| Silica-based Sorbent | Reversed phase | Normal Phase | Weak Ion Exchange | Strong Ion Exchange |
|-----------------------|--|---|---|---------------------|
| | C18 C18-U C8 C2 PH | Silica CN (can also be used as reversed phase) | NH ₂ PSA | SCX SAX |
| Polymer-based Sorbent | Universal Phase | Mix Mode* | Mix Mode** | |
| | PLS | PXC PXA | PWC PWA | |
| Other | Inorganic Material | | Double Layer Sorbent | |
| | Al ₂ O ₃ (Acidic) Al ₂ O ₃ (Neutral) Al ₂ O ₃ (Basic) Florisil (Magnesium silicate) CARB (Graphitized carbon black) Na ₂ SO ₄ | | CARB / NH ₂ CARB / PSA SAX / PSA | |

*Reversed phase and strong ion exchange

**Reversed phase and weak ion exchange

Brand Cross Reference

| Dikma | Waters | Agilent | Supelco |
|--------------------------|-------------------------|--------------------------|--------------------|
| ProElut™ PLS | Oasis HLB | Plexa | - |
| ProElut™ PXC | Oasis MCX | Plexa PCX | - |
| ProElut™ PXA | Oasis MAX | - | - |
| ProElut™ PWC | Oasis WCX | - | - |
| ProElut™ PWA | Oasis WAX | - | - |
| ProElut™ C18 | Sep-pak C18 | BondElut C18 | ENVI-18, LC-18 |
| ProElut™ C18-U | - | BondElut C18-OH | - |
| ProElut™ C8 | Sep-pak C8 | BondElut C8 | - |
| ProElut™ C2 | - | BondElut C2 | - |
| ProElut™ PH | - | BondElut PH | LC-PH |
| ProElut™ CN | Sep-pak CN | BondElut CN | LC-CN |
| ProElut™ NH ₂ | Sep-pak NH ₂ | BondElut NH ₂ | LC-NH ₂ |
| ProElut™ PSA | - | BondElut PSA | - |
| ProElut™ Silica | Sep-pak Silica | BondElut Silica | LC-Silica |
| ProElut™ SCX | - | BondElut SCX | LC-SCX |
| ProElut™ SAX | - | BondElut SAX | LC-SAX |

Technical Reference

A. Typical elution sequence of PLS, PXC, PXA, PWA and PWC

ProElut™ PLS (60 mg / 3 mL)

| Suggested Method | Purpose |
|--|---|
| Condition / equilibrate 3 mL MeOH / 3 mL H ₂ O | Ready for use |
| Load sample solution 3 mL | The compounds of interest will be adsorbed by the sorbent |
| Wash 3 mL H ₂ O (5% MeOH) | To remove aqueous soluble materials and disruptors |
| Elution 3 mL MeOH | To get compounds of interest that previously adsorbed on the sorbent by non-polar interaction |
| Evaporate and reconstitute | For HPLC or GC analysis |

ProElut™ PXC (60 mg / 3 mL)

| Suggested Method | Purpose |
|--|--|
| Condition / equilibrate 3 mL MeOH / 3 mL H ₂ O | Ready for use |
| Load acidified sample solution 3 mL | Protonated basic compounds (under low pH) will approach sulfonic group by Coulomb force The neutral and acidic compounds will be adsorbed on the sorbent by non-polar interaction |
| Wash 1 3 mL 0.1 M HCl | To remove aqueous soluble materials and disruptors |
| Wash 2 3 mL MeOH | To remove compounds that adsorbed on the sorbent by non-polar interaction |
| Elution 3 mL MeOH (5% NH ₄ OH) | Neutralize the basic compounds that adsorbed on the sorbent by Coulomb force and carry them out |
| Evaporate and reconstitute | For HPLC or GC analysis |

ProElut™ PXA (60 mg / 3 mL)

| Suggested Method | Purpose |
|--|--|
| Condition / equilibrate 3 mL MeOH / 3 mL H ₂ O | Ready for use |
| Load sample solution 3 mL | Negative charged acidic compounds will approach quaternary amino group by Coulomb force The neutral and acidic compounds will be adsorbed on the sorbent by non-polar interaction |
| Wash 1 3 mL 5% NH ₄ OH | To remove aqueous soluble materials and disruptors, including salts and proteins The interaction between acid compounds and the quaternary amino group is reinforced |
| Wash 2 3 mL MeOH | To remove compounds that adsorbed on the sorbent by non-polar interaction |
| Elution 3 mL MeOH (2% HCOOH) | Neutralize the acidic compounds that adsorbed on the sorbent by Coulomb force and carry them out |
| Evaporate and reconstitute | For HPLC or GC analysis |

ProElut™ SPE

ProElut™ PWC (60 mg / 3 mL)

| Suggested Method | Purpose |
|---|---|
| Condition / equilibrate 3 mL MeOH / 3 mL 5% NH ₄ OH | Ready for use Add ammonia to make the carboxyl functional group negatively charged |
| Load sample solution 3 mL | Protonated strong basic compounds will approach carboxyl group by Coulomb force The neutral and weak- / mid- basic compounds will be absorbed on the sorbent by non-polar interaction |
| Wash 1 3 mL 5% NH ₄ OH | To remove aqueous soluble materials and disruptors, including salts and proteins The interaction between strong basic compounds and the carboxyl group is reinforced |
| Wash 2 3 mL MeOH | To remove compounds that adsorbed on the sorbent by non-polar interaction |
| Elution 3 mL MeOH (2% HCOOH) | To neutralize the carboxyl (negatively charged) so that the Coulomb force between analyte and functional group is cut off and therefore the strong basic compounds will be carried away by methanol |
| Evaporate and reconstitute | For HPLC or GC analysis |

ProElut™ PWA (60 mg / 3 mL)

| Suggested Method | Purpose |
|--|---|
| Condition / equilibrate 3 mL MeOH / 3 mL 2% HCOOH | Ready for use Add formic acid to make the piperazine functional group protonated |
| Load sample solution 3 mL | Negatively charged strong acidic compounds will approach sorbent functional group by Coulomb force The neutral and weak- / mid- basic compounds will be absorbed on the sorbent by non-polar interaction |
| Wash 1 3 mL 2% HCOOH | To remove aqueous soluble materials and disruptors, including salts and proteins The interaction between strong acidic compounds and the piperazine group is reinforced |
| Wash 2 3 mL MeOH | To remove compounds that adsorbed on the sorbent by non-polar interaction |
| Elution 3 mL MeOH (5% NH ₄ OH) | To neutralize the protonated functional group so that the Coulomb force between analyte and functional group is cut off and therefore the strong acidic compounds will be carried away by methanol |
| Evaporate and reconstitute | For HPLC or GC analysis |

B. ProElut™ SPE Sorbents Weight Based on Sample Size

| Bed Weight (mg) | Bed Capacity (mg)* | Minimum Elution Volume (μL) |
|-----------------|--------------------|-----------------------------|
| 50 | 2.5 | 125 |
| 100 | 5 | 250 |
| 150 | 7.5 | 375 |
| 200 | 10 | 500 |
| 500 | 25 | 1250 |
| 1000 | 50 | 2500 |

*This value depends on the analyte and sample matrix. As a rule of thumb, the bed capacity can be estimated with ~5% of the bed weight.

ProElut™ SPE Sorbents

Silica-based Sorbents

ProElut™ C18

ProElut™ C18 is the most broadly used SPE cartridge. It can be used to adsorb non-polar, slightly polar and mid-polar compounds. Polar materials such as salt cannot be retained on the sorbent, which makes ProElut™ C18 an excellent choice for desalting samples. In addition, non-polar and slightly polar disruptors in matrix such as fats, PAHs and phthalates can be retained by the sorbent, leaving ionic analytes eluted in the collector for reconstitution.

| | |
|---------------------|--|
| Base material | Spherical silica, Particle Size: 50 μm, Pore Size: 60 Å Specific Surface Area: 500 m ² /g |
| Functional group | $\begin{array}{c} \text{CH}_3 \\ \\ \text{---Si---C}_{18}\text{H}_{37} \\ \\ \text{CH}_3 \end{array}$ |
| Endcapping | Yes |
| Carbon loading | 19% |
| Retention mechanism | Reversed phase |
| Application | For reversed phase extraction of non-polar to moderately polar compounds, such as antibiotics, barbiturates, benzodiazepines, caffeine, drugs, dyes, essential oils, fat-soluble vitamins, fungicides, herbicides, pesticides, hydrocarbons, parabens, phenols, phthalate esters, steroids, surfactants, theophylline and water-soluble vitamins |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|-----------------------------------|--------|-----|----------|---------|
| ProElut™ C18 SPE Tubes | | | | |
| 50 mg | 1 mL | 100 | 63101 | \$125 |
| 100 mg | 1 mL | 100 | 63102 | 125 |
| 200 mg | 3 mL | 50 | 63103 | 95 |
| 500 mg | 3 mL | 50 | 63104 | 115 |
| 500 mg | 6 mL | 30 | 63105 | 80 |
| 1 g | 6 mL | 30 | 63106 | 105 |
| 2 g | 12 mL | 20 | 63107 | 120 |
| 5 g | 20 mL | 20 | 63108 | 190 |
| 10 g | 60 mL | 10 | 63109 | 260 |
| ProElut™ C18 Bulk Sorbents | | | | |
| 10 g | | 1 | 63181 | Inquire |
| 100 g | | 1 | 63182 | Inquire |
| 1 kg | | 1 | 63183 | Inquire |

ProElut™ C18-U

ProElut™ C18-U has a non-encapped octadecyl bonded phase that enables the silanols on the silica surface to be more active. However, the well-controlled silanol activity exhibits excellent selectivity towards polar compounds, especially amines such as tetracyclines.

| | |
|---------------------|--|
| Base material | Spherical silica, Particle Size: 50 μm, Pore Size: 60 Å, Specific Surface Area: 500 m ² /g |
| Functional group | $\begin{array}{c} \text{CH}_3 \\ \\ \text{---Si---C}_{18}\text{H}_{37} \\ \\ \text{CH}_3 \\ \\ \text{---Si---OH} \\ \end{array}$ |
| Endcapping | No |
| Carbon loading | 17% |
| Retention mechanism | Reversed phase |
| Application | Similar to ProElut™ C18, but with enhanced retention for polar compounds |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|-------------------------------------|--------|-----|----------|---------|
| ProElut™ C18-U SPE Tubes | | | | |
| 50 mg | 1 mL | 100 | 63501 | \$125 |
| 100 mg | 1 mL | 100 | 63502 | 125 |
| 200 mg | 3 mL | 50 | 63503 | 95 |
| 500 mg | 3 mL | 50 | 63504 | 115 |
| 500 mg | 6 mL | 30 | 63505 | 80 |
| 1 g | 6 mL | 30 | 63506 | 105 |
| 2 g | 12 mL | 20 | 63507 | 120 |
| 5 g | 20 mL | 20 | 63508 | 190 |
| 10 g | 60 mL | 10 | 63509 | 260 |
| ProElut™ C18-U Bulk Sorbents | | | | |
| 10 g | | 1 | 63581 | Inquire |
| 100 g | | 1 | 63582 | Inquire |
| 1 kg | | 1 | 63583 | Inquire |

ProElut™ SPE

ProElut™ C8

ProElut™ C8 is very similar to the C18 phase, but has a shorter chain. This makes the phase less non-polar than C18, leaving non-polar compounds less retained by the sorbents. In this case, those compounds retained too strongly on the C18 can be effectively eluted if you choose the C8 phase. In addition, the C8 phase for polar interaction is somewhat higher than C18 because there is less coverage of the silica surface. However, this polar interaction is not the main characteristic of C8 phase.

| | |
|---------------------|---|
| Base material | Spherical silica, Particle Size: 50 μm, Pore Size: 60 Å Specific Surface Area: 500 m ² /g |
| Functional group | $\begin{array}{c} \text{CH}_3 \\ \\ \text{---Si---C}_8\text{H}_{17} \\ \\ \text{CH}_3 \end{array}$ |
| Endcapping | Yes |
| Carbon loading | 11% |
| Retention mechanism | Reversed phase |
| Application | For reversed phase extraction of non-polar to moderately polar compounds, such as antibiotics, barbiturates, benzodiazepines, caffeine, drugs, dyes, essential oils, fat-soluble vitamins, fungicides, herbicides, pesticides, hydrocarbons, parabens, phenols, phthalate esters, steroids, surfactants, theophylline, and water-soluble vitamins |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|----------------------------------|--------|-----|----------|---------|
| ProElut™ C8 SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 63702 | \$125 |
| 200 mg | 3 mL | 50 | 63703 | 95 |
| 500 mg | 3 mL | 50 | 63704 | 115 |
| 500 mg | 6 mL | 30 | 63705 | 80 |
| 1 g | 6 mL | 30 | 63706 | 105 |
| 2 g | 12 mL | 20 | 63707 | 120 |
| 5 g | 20 mL | 20 | 63708 | 190 |
| 10 g | 60 mL | 10 | 63709 | 260 |
| ProElut™ C8 Bulk Sorbents | | | | |
| 10 g | | 1 | 63781 | Inquire |
| 100 g | | 1 | 63782 | Inquire |
| 1 kg | | 1 | 63783 | Inquire |

ProElut™ C2

ProElut™ C2 has the shortest carbon chain among non-polar silica-based phases. It provides the least retentive ability for non-polar compounds and somewhat higher polar interaction than C8 and C18 phases. The polarity of C2 is slightly lower than a cyano phase for polar interactions.


| | |
|---------------------|---|
| Base material | Spherical silica, Particle Size: 50 μm, Pore Size: 60 Å Specific Surface Area: 500 m ² /g |
| Functional group | $\begin{array}{c} \text{CH}_3 \\ \\ \text{---Si---C}_2\text{H}_5 \\ \\ \text{CH}_3 \end{array}$ |
| Endcapping | Yes |
| Carbon loading | 5.6% |
| Retention mechanism | Reversed phase |
| Application | Plasma, urine, aqueous samples |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|----------------------------------|--------|-----|----------|---------|
| ProElut™ C2 SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 65602 | \$125 |
| 500 mg | 3 mL | 50 | 65604 | 115 |
| 1 g | 6 mL | 30 | 65606 | 105 |
| 2 g | 12 mL | 20 | 65607 | 120 |
| 5 g | 20 mL | 20 | 65608 | 190 |
| 10 g | 60 mL | 10 | 65609 | 260 |
| ProElut™ C2 Bulk Sorbents | | | | |
| 10 g | | 1 | 65681 | Inquire |
| 100 g | | 1 | 65682 | Inquire |
| 1 kg | | 1 | 65683 | Inquire |

ProElut™ PH

ProElut™ PH has a similar polarity as that of C8 phase. However, it shows different selectivity to other non-polar phases due to the presence of conjugated double bonds. In addition, the planar shape of benzene and its electron distribution make it much more retentive to aromatic compounds.


| | |
|---------------------|--|
| Base material | Spherical silica, Particle Size: 50 μm , Pore Size: 60 \AA Specific Surface Area: 500 m^2/g |
| Functional group |  |
| Endcapping | Yes |
| Carbon loading | 10% |
| Retention mechanism | Reversed phase |
| Application | Volatiles in water: PAHs, PAEs, PCBs, pesticides, herbicides, phenols Biological fluids: blood, urine |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|----------------------------------|--------|-----|----------|---------|
| ProElut™ PH SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 63902 | \$125 |
| 500 mg | 3 mL | 50 | 63904 | 115 |
| 1 g | 6 mL | 30 | 63906 | 105 |
| 2 g | 12 mL | 20 | 63907 | 120 |
| 5 g | 20 mL | 20 | 63908 | 190 |
| 10 g | 60 mL | 10 | 63909 | 260 |
| ProElut™ PH Bulk Sorbents | | | | |
| 10 g | | 1 | 63981 | Inquire |
| 100 g | | 1 | 63982 | Inquire |
| 1 kg | | 1 | 63983 | Inquire |

ProElut™ CN

ProElut™ CN is a mid-polar phase SPE column. We recommend using it to extract samples for which analytes are irreversibly retained on C8 and C18 phases.


| | |
|---------------------|--|
| Base material | Spherical silica, Particle Size: 50 μm , Pore Size: 60 \AA Specific Surface Area: 500 m^2/g |
| Functional group |  |
| Endcapping | Yes |
| Carbon loading | 8% |
| Retention mechanism | Reversed phase or normal phase |
| Application | Pesticides in water, metabolites |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|----------------------------------|--------|-----|----------|---------|
| ProElut™ CN SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 63802 | \$125 |
| 500 mg | 3 mL | 50 | 63804 | 115 |
| 1 g | 6 mL | 30 | 63806 | 105 |
| 2 g | 12 mL | 20 | 63807 | 120 |
| 5 g | 20 mL | 20 | 63808 | 190 |
| 10 g | 60 mL | 10 | 63809 | 260 |
| ProElut™ CN Bulk Sorbents | | | | |
| 10 g | | 1 | 63881 | Inquire |
| 100 g | | 1 | 63882 | Inquire |
| 1 kg | | 1 | 63883 | Inquire |

ProElut™ Silica

ProElut™ Silica is the most polar SPE sorbent. It is very effective for separating compounds with similar structures and extracting polar compounds in non-polar solvents. In addition, the silica surface silanols have slight anion exchange properties that can be used to remove organic acids and phenols in extracts.

| | |
|---------------------|---|
| Base material | Spherical silica, Particle Size: 50 μm , Pore Size: 60 \AA Specific Surface Area: 500 m^2/g |
| Functional group |  |
| Endcapping | No |
| Carbon loading | - |
| Retention mechanism | Normal phase or weak anion exchange |
| Application | Extract polar compounds from non-polar matrix Remove polar hydrocarbons, organic acids, and phenols in extracts. Separate compounds with very similar structures (isomers) |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|--------------------------------------|--------|-----|----------|---------|
| ProElut™ Silica SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 63002 | \$125 |
| 500 mg | 3 mL | 50 | 63004 | 115 |
| 1 g | 6 mL | 30 | 63006 | 105 |
| ProElut™ Silica Bulk Sorbents | | | | |
| 10 g | | 1 | 63081 | Inquire |
| 100 g | | 1 | 63082 | Inquire |
| 1 kg | | 1 | 63083 | Inquire |

ProElut™ SPE

ProElut™ NH₂

ProElut™ NH₂ has both polar and weak anion exchange interactions. It can effectively absorb compounds with a polar functional group (-OH, -NH₂, -SH, etc.) by hydrogen bonding from non-polar solvents such as hexane. In addition, it has weaker anion exchange property than SAX (a quaternary amine sorbent that is always charged) and is therefore an excellent choice for retention of very strong anions that are always irreversibly adsorbed on a SAX sorbent, such as sulfonic acid.

| | |
|---------------------|---|
| Base material | Spherical silica, Particle Size: 50 μm, Pore Size: 60 Å Specific Surface Area: 500 m ² /g |
| Functional group | $\begin{array}{c} \text{CH}_3 \\ \\ \text{---Si---}(\text{CH}_2)_3\text{NH}_2 \\ \\ \text{CH}_3 \end{array}$ |
| Endcapping | No |
| Carbon loading | 5.5% |
| Retention mechanism | Normal phase or anion exchange |
| Application | Extract polar compounds from non-polar matrix Remove polar hydrocarbons, organic acids, and phenols in extracts Separate compounds with very similar structures (isomers) |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|--|--------|-----|----------|---------|
| ProElut™ NH₂ SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 63302 | \$125 |
| 200 mg | 3 mL | 50 | 63303 | 95 |
| 500 mg | 3 mL | 50 | 63304 | 115 |
| 500 mg | 6 mL | 30 | 63305 | 80 |
| 1 g | 6 mL | 30 | 63306 | 105 |
| 2 g | 12 mL | 20 | 63307 | 120 |
| 5 g | 20 mL | 20 | 63308 | 190 |
| 10 g | 60 mL | 10 | 63361 | 260 |
| ProElut™ NH₂ Bulk Sorbents | | | | |
| 10 g | | 1 | 63381 | Inquire |
| 100 g | | 1 | 63382 | Inquire |
| 1 kg | | 1 | 63383 | Inquire |

ProElut™ PSA

ProElut™ PSA sorbent contains two different amino groups, one primary and one secondary. It gives comparatively higher pK_a and ionic capacity relative to ProElut™ NH₂. The PSA sorbent is an excellent choice for extracting polar compounds from non-polar solvents. The compounds that are retained too strongly on a NH₂ sorbent can be effectively eluted on a PSA sorbent. In addition, the PSA functional group is a very effective bidentate ligand in chelation applications.

| | |
|---------------------|---|
| Base material | Spherical silica, Particle Size: 50 μm, Pore Size: 60 Å Specific Surface Area: 500 m ² /g |
| Functional group | $\begin{array}{c} \text{CH}_3 \\ \\ \text{---Si---}(\text{CH}_2)_3\text{NH}(\text{CH}_2)_2\text{NH}_2 \\ \\ \text{CH}_3 \end{array}$ |
| Endcapping | No |
| Carbon loading | 8.5% |
| Retention mechanism | Normal phase or weak anion exchange |
| Application | Extract polar compounds from non-polar matrix Remove polar hydrocarbons, organic acids, and phenols in extracts Separate compounds with very similar structures (isomers) |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|-----------------------------------|--------|-----|----------|---------|
| ProElut™ PSA SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 63202 | \$150 |
| 200 mg | 3 mL | 50 | 63203 | 115 |
| 500 mg | 3 mL | 50 | 63204 | 140 |
| 500 mg | 6 mL | 30 | 63205 | 115 |
| 1 g | 6 mL | 30 | 63206 | 160 |
| 2 g | 12 mL | 20 | 63207 | 188 |
| 5 g | 20 mL | 20 | 63208 | 223 |
| 10 g | 60 mL | 10 | 63209 | 307 |
| ProElut™ PSA Bulk Sorbents | | | | |
| 10 g | | 1 | 63281 | Inquire |
| 100 g | | 1 | 63282 | Inquire |
| 1 kg | | 1 | 63283 | Inquire |

ProElut™ SCX

ProElut™ SCX sorbent has benzenesulfonic acid as a bonded functional group with a very low pK_a . The presence of the benzene ring in the functional group increases its potential for non-polar interaction. The two properties are quite useful in the absorption of cationic organic compounds from aqueous systems where non-polar compounds are seen.

| | |
|---------------------|--|
| Base material | Spherical silica, Particle Size: 50 μm , Pore Size: 60 \AA Specific Surface Area: 500 m^2/g |
| Functional group | $\begin{array}{c} \text{CH}_3 \\ \\ \text{---Si---}(\text{CH}_2)_3\text{C}_6\text{H}_4\text{SO}_3^-\text{H}^+ \\ \\ \text{CH}_3 \end{array}$ |
| Endcapping | No |
| Carbon loading | 10.9% |
| Retention mechanism | Strong cation exchange |
| Application | Basic compounds in aqueous solution |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|-----------------------------------|--------|-----|----------|---------|
| ProElut™ SCX SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 63602 | \$150 |
| 500 mg | 3 mL | 50 | 63604 | 140 |
| 500 mg | 6 mL | 30 | 63606 | 115 |
| 2 g | 12 mL | 20 | 63607 | 147 |
| 5 g | 20 mL | 20 | 63608 | 211 |
| 10 g | 60 mL | 10 | 63609 | 271 |
| ProElut™ SCX Bulk Sorbents | | | | |
| 10 g | | 1 | 63681 | Inquire |
| 100 g | | 1 | 63682 | Inquire |
| 1 kg | | 1 | 63683 | Inquire |

ProElut™ SAX

ProElut™ SAX sorbent has trimethylaminopropyl as a bonded functional group with a very high pK_a . The presence of the benzene ring in the functional group increases its potential for non-polar interaction. The two properties are quite useful in the absorption of anionic organic compounds from aqueous systems where non-polar compounds are seen.

| | |
|---------------------|---|
| Base material | Spherical silica, Particle Size: 50 μm , Pore Size: 60 \AA Specific Surface Area: 500 m^2/g |
| Functional group | $\begin{array}{c} \text{CH}_3 \\ \\ \text{---Si---}(\text{CH}_2)_3\text{N}^+(\text{CH}_3)_3\text{Cl}^- \\ \\ \text{CH}_3 \end{array}$ |
| Endcapping | No |
| Carbon loading | 8% |
| Retention mechanism | Strong anion exchange |
| Application | Carboxylic acids in aqueous solution |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|-----------------------------------|--------|-----|----------|---------|
| ProElut™ SAX SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 63402 | \$150 |
| 500 mg | 3 mL | 50 | 63404 | 140 |
| 500 mg | 6 mL | 30 | 63406 | 115 |
| 2 g | 12 mL | 20 | 63407 | 147 |
| 5 g | 20 mL | 20 | 63408 | 211 |
| 10 g | 60 mL | 10 | 63409 | 271 |
| ProElut™ SAX Bulk Sorbents | | | | |
| 10 g | | 1 | 63481 | Inquire |
| 100 g | | 1 | 63482 | Inquire |
| 1 kg | | 1 | 63483 | Inquire |

ProElut™ SPE

Polymer-based Sorbents

ProElut™ PLS-Hydrophilic-Lipophilic-Balance Copolymer, Reversed Phase Sorbent

ProElut™ PLS is a hydrophilic polystyrene / divinylbenzene copolymer sorbent, designed to expand the SPE application fields and improve extraction efficiency. This sorbent contains the lipophilic divinylbenzene and the hydrophilic pyrrolidone. The hydrophilic-lipophilic balance is a reversed phase sorbent maintaining retention for non-polar and polar analytes. Compared to traditional silica-based reversed phase sorbent (C18), ProElut™ PLS features are as follows:

(1) Real Versatility

- High retention of hydrophilic compounds and lipophilic compounds
- Applications covering the non-polar, weakly polar, and polar compounds to overcome the C18 sorbent poor retention of polar compounds

(2) Higher Stability

- PLS is water-wettable
- Maintaining high retention and capacity after activation even if the SPE cartridge runs dry

(3) Wider pH Range

PLS is a polymer sorbent that can be used in the range of pH 0 - 14 while the silica-based sorbents can only be used in the range of pH 2 - 7.5.


(4) Higher Capacity

The PLS sorbent has greater capacity for more compounds. It reduces breakthrough potential and improves reproducibility.

(5) No Secondary Interaction


The residue silanols of silica-based sorbents can adsorb basic analytes resulting in low recovery. PLS is a polymer sorbent and there is no silanol activity leading to high recovery.

Structure of ProElut™ Polymer-based Sorbents



ProElut™ PLS

ProElut™ PLS is a highly cross-linked polystyrene-divinylbenzene (PS-DVB) copolymer with high surface area (800 m²/g) and high capacity. It is an excellent choice for extraction of polar analytes in aqueous solvents where traditional C18 and C8 sorbents are not advisable because they are not "wetttable". It is ideal in screening applications where a broad range of analytes can be extracted.


| | |
|----------------------------------|--|
| Base material | Porous, highly cross-linked, spherical PS-DVB, Particle Size: 50 μm, Pore Size: 80 Å, Specific Surface Area: 800 m ² /g |
| Functional group (Pyrrolidinone) |  |
| Retention mechanism | Non-polar and polar interactions |
| Application | Pharmaceutical residues in animal tissue, such as tetracyclines, chloromycetin, sulfonamides, abamectin, macrolide antibiotics, nitrofurans, and pesticides in vegetables Environmental samples, such as PAHs, PAEs, phenols, and endocrine disruptors Biological samples, pharmaceuticals, and metabolites in plasma, serum, or urine |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|-------------------------------|--------|-----|----------|-------|
| ProElut™ PLS SPE Tubes | | | | |
| 30 mg | 1 mL | 100 | 68002 | \$159 |
| 60 mg | 1 mL | 100 | 68011 | 184 |
| 60 mg | 3 mL | 50 | 68003 | 221 |
| 150 mg | 6 mL | 30 | 68004 | 136 |
| 200 mg | 6 mL | 30 | 68012 | 155 |
| 500 mg | 6 mL | 30 | 68005 | 195 |
| 500 mg | 12 mL | 20 | 68007 | 156 |
| 1 g | 20 mL | 20 | 68008 | 242 |
| 6 g | 60 mL | 10 | 68009 | 491 |

ProElut™ PXC

ProElut™ PXC is a highly cross-linked PS-DVB copolymer with sulfonic acid as the functional group. It has both non-polar and cation exchange interactions and is therefore an excellent choice for extraction of basic organic compounds.


| | |
|----------------------------------|--|
| Base material | Porous, highly cross-linked, spherical PS-DVB, Particle Size: 50 μm, Pore Size: 80 Å, Specific Surface Area: 800 m ² /g |
| Functional group (Sulfonic acid) |  |
| Retention mechanism | Non-polar interaction and cation exchange |
| Application | Basic compounds, such as sulfonamides and clenbuterol Biological samples, pharmaceuticals, and metabolites in plasma, serum, or urine |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|-------------------------------|--------|-----|----------|-------|
| ProElut™ PXC SPE Tubes | | | | |
| 30 mg | 1 mL | 100 | 68202 | \$143 |
| 60 mg | 3 mL | 50 | 68203 | 198 |
| 150 mg | 6 mL | 30 | 68204 | 122 |
| 200 mg | 6 mL | 30 | 68212 | 139 |
| 500 mg | 12 mL | 20 | 68207 | 140 |
| 1 g | 20 mL | 20 | 68208 | 218 |
| 6 g | 60 mL | 10 | 68209 | 442 |

ProElut™ PXA

ProElut™ PXA is a highly cross-linked PS-DVB copolymer with a quaternary amino as the functional group. It has both non-polar and anion exchange interactions and is therefore an excellent choice for extraction of acidic organic compounds, especially those containing carboxyl and phenolic hydroxyl.

| | |
|-------------------------------------|--|
| Base material | Porous, highly cross-linked, spherical PS-DVB, Particle Size: 50 μm, Pore Size: 80 Å, Specific Surface Area: 800 m ² /g |
| Functional group (Quaternary amino) |  |
| Retention mechanism | Non-polar interaction and anion exchange |
| Application | Compounds with groups as carboxyl and phenolic hydroxyl |


Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|-------------------------------|--------|-----|----------|-------|
| ProElut™ PXA SPE Tubes | | | | |
| 30 mg | 1 mL | 100 | 68302 | \$159 |
| 60 mg | 3 mL | 50 | 68303 | 221 |
| 150 mg | 6 mL | 30 | 68304 | 136 |
| 500 mg | 12 mL | 20 | 68307 | 156 |
| 1 g | 20 mL | 20 | 68308 | 242 |
| 6 g | 60 mL | 10 | 68309 | 491 |

ProElut™ SPE

ProElut™ PWC

ProElut™ PWC is a highly cross-linked PS-DVB copolymer with carboxyl as the functional group. It has both non-polar and weak cation exchange interactions and is therefore an excellent choice for extraction of strong basic compounds.


| | |
|-----------------------------|--|
| Base material | Porous, highly cross-linked, spherical PS-DVB, Particle Size: 50 μm , Pore Size: 80 \AA , Specific Surface Area: 800 m^2/g |
| Functional group (Carboxyl) |  |
| Retention mechanism | Non-polar interaction and weak cation exchange |
| Application | Strong basic compounds, quaternary ammonium salts |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|-------------------------------|--------|-----|----------|-------|
| ProElut™ PWC SPE Tubes | | | | |
| 30 mg | 1 mL | 100 | 65711 | \$159 |
| 60 mg | 3 mL | 50 | 65712 | 221 |
| 150 mg | 6 mL | 30 | 65713 | 136 |

ProElut™ PWA

ProElut™ PWA is a highly cross-linked PS-DVB copolymer with piperazine as the functional group. It has both non-polar and weak anion exchange interactions and is therefore an excellent choice for extraction of strong acidic compounds.

| | |
|-------------------------------|--|
| Base material | Porous, highly cross-linked, spherical PS-DVB, Particle Size: 50 μm , Pore Size: 80 \AA , Specific Surface Area: 800 m^2/g |
| Functional group (Piperazine) |  |
| Retention mechanism | Non-polar interaction and weak anion exchange |
| Application | For purification of strong acidic compounds |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|-------------------------------|--------|-----|----------|-------|
| ProElut™ PWA SPE Tubes | | | | |
| 30 mg | 1 mL | 100 | 65811 | \$159 |
| 60 mg | 3 mL | 50 | 65812 | 221 |
| 150 mg | 6 mL | 30 | 65813 | 136 |

Specific Sorbents

ProElut™ Florisil

Florisil is a highly selective adsorbent that has extensive utility in sample preparation, preparative and analytical chromatography. This sorbent is unique because it is comprised of extremely white, hard-powdered synthetic magnesium-silica gel.

| | |
|---------------------|---|
| Base material | Magnesium silicate, Particle Size: 150 - 200 μm |
| Functional group | MgSiO ₃ |
| Retention mechanism | Polar interaction |
| Application | For extraction of drugs, dyes, herbicides, pesticides, nitrogen compounds, organic acids, phenols, steroids, PCBs, and PAHs |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|--|--------|-----|----------|---------|
| ProElut™ Florisil SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 65002 | \$125 |
| 500 mg | 3 mL | 50 | 65004 | 115 |
| 500 mg | 6 mL | 30 | 65005 | 80 |
| 1 g | 6 mL | 30 | 65006 | 105 |
| 2 g | 12 mL | 20 | 65007 | 120 |
| 5 g | 20 mL | 20 | 65008 | 190 |
| 10 g | 60 mL | 10 | 65009 | 260 |
| ProElut™ Florisil Bulk Sorbents | | | | |
| 10 g | | 1 | 65081 | Inquire |
| 100 g | | 1 | 65082 | Inquire |
| 1 kg | | 1 | 65083 | Inquire |

ProElut™ AL-A

| | |
|---------------------|---|
| Base material | Al ₂ O ₃ (Acidic), Particle Size: 125 μm |
| pH | 4.5 |
| Retention mechanism | Lewis acid, polar, and ion exchange interactions |
| Application | Polar and anionic compounds |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|------------------------------------|--------|-----|----------|---------|
| ProElut™ AL-A SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 65102 | \$125 |
| 500 mg | 3 mL | 50 | 65104 | 115 |
| 1 g | 6 mL | 30 | 65106 | 105 |
| 2 g | 12 mL | 20 | 65107 | 120 |
| 5 g | 20 mL | 20 | 65108 | 190 |
| 10 g | 60 mL | 10 | 65109 | 260 |
| ProElut™ AL-A Bulk Sorbents | | | | |
| 10 g | | 1 | 65181 | Inquire |
| 100 g | | 1 | 65182 | Inquire |
| 1 kg | | 1 | 65183 | Inquire |

ProElut™ CARB

| | |
|---------------|---|
| Base material | Graphitized carbon black, Particle Size: 120 - 400 μm |
| Application | In agriculture residues analysis, used to remove pigments in fruits and vegetables For purification of samples such as groundwater, fruits, and vegetables |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|--------------------------------|--------|-----|----------|-------|
| ProElut™ CARB SPE Tubes | | | | |
| 250 mg | 3 mL | 50 | 65403 | \$200 |
| 500 mg | 6 mL | 30 | 65405 | 210 |

ProElut™ AL-N

| | |
|---------------------|---|
| Base material | Al ₂ O ₃ (Neutral), Particle Size: 125 μm |
| pH | 7.5 |
| Retention mechanism | Lewis acid, polar, and ion exchange interactions |
| Application | Extract polar compounds from non-polar matrix Remove polar hydrocarbons, organic acids, and phenols in extracts. Separate compounds with very similar structures (isomers) |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|------------------------------------|--------|-----|----------|---------|
| ProElut™ AL-N SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 65302 | \$125 |
| 500 mg | 3 mL | 50 | 65304 | 115 |
| 1 g | 6 mL | 30 | 65306 | 105 |
| 2 g | 12 mL | 20 | 65307 | 120 |
| 5 g | 20 mL | 20 | 65308 | 190 |
| 10 g | 60 mL | 10 | 65309 | 260 |
| ProElut™ AL-N Bulk Sorbents | | | | |
| 10 g | | 1 | 65381 | Inquire |
| 100 g | | 1 | 65382 | Inquire |
| 1 kg | | 1 | 65383 | Inquire |

ProElut™ AL-B

| | |
|---------------------|--|
| Base material | Al ₂ O ₃ (Basic), Particle Size: 125 μm |
| pH | 10.0 |
| Retention mechanism | Lewis acid, polar, and ion exchange interactions |
| Application | Polar and anionic compounds, amines |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|------------------------------------|--------|-----|----------|---------|
| ProElut™ AL-B SPE Tubes | | | | |
| 100 mg | 1 mL | 100 | 65202 | \$125 |
| 500 mg | 3 mL | 50 | 65204 | 115 |
| 1 g | 6 mL | 30 | 65206 | 105 |
| 2 g | 12 mL | 20 | 65207 | 120 |
| 5 g | 20 mL | 20 | 65208 | 190 |
| 10 g | 60 mL | 10 | 65209 | 260 |
| ProElut™ AL-B Bulk Sorbents | | | | |
| 10 g | | 1 | 65281 | Inquire |
| 100 g | | 1 | 65282 | Inquire |
| 1 kg | | 1 | 65283 | Inquire |

ProElut™ CARB / NH₂

| | |
|---------------|---|
| Base material | Packing of equivalent volume of CARB and amine (NH ₂) |
| Application | Widely used in analysis of pesticide residues (many different varieties) in foods To remove pigments, fatty acids, and phenols from analytes |

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|---|--------|-----|----------|-------|
| ProElut™ CARB / NH₂ SPE Tubes | | | | |
| 500 mg / 500 mg | 6 mL | 30 | 64105 | \$240 |

ProElut™ SPE

ProElut™ GLASS SPE Tube

ProElut™ glass cartridges are designed for high-purity extraction as the inert glass body completely eliminates the pollution from plasticizers, such as phthalates. ProElut™ glass SPE tubes are standard size, with high quality sorbent and special purification frits, to assure stability and reproducibility.

Ordering Information

| Mass | Volume | Qty | Cat. No. | Price |
|--|--------|-----|----------|-------|
| ProElut™ PLS GLASS SPE Tubes | | | | |
| 200 mg | 6 mL | 30 | 68012G | \$189 |
| 500 mg | 6 mL | 30 | 68005G | 397 |
| ProElut™ C18 GLASS SPE Tubes | | | | |
| 500 mg | 6 mL | 30 | 63105G | 228 |
| 1 g | 6 mL | 30 | 63106G | 302 |
| ProElut™ Florisil GLASS SPE Tubes | | | | |
| 500 mg | 6 mL | 30 | 65005G | 271 |
| 1 g | 6 mL | 30 | 65006G | 302 |
| ProElut™ PSA GLASS SPE Tubes | | | | |
| 500 mg | 6 mL | 30 | 63205G | 200 |
| 1 g | 6 mL | 30 | 63206G | 270 |
| ProElut™ Silica GLASS SPE Tubes | | | | |
| 500 mg | 6 mL | 30 | 63005G | 198 |
| 1 g | 6 mL | 30 | 63006G | 254 |

