

Immobilized metal ion affinity chromatography


Immobilized metal ion affinity chromatography (IMAC) separates most proteins with exposed histidine, cysteine and tryptophan on their surface. IMAC is an excellent technique for optimization and purification of His-tagged proteins. The technique is ideal for capture directly from clarified cell lysate. The target protein is collected in a highly purified and concentrated form.

Several factors influence the final purity of a His-tagged protein after an IMAC purification, for example, position of the tag (C- or N-terminal), the length of the tag, immobilized metal ion (Ni^{2+} , Co^{2+} , Zn^{2+} , Cu^{2+}) and the ligand immobilized on the matrix (NTA or IDA). To make the optimization of His-tagged protein purifications as efficient as possible, Bio-Works offers products with many combinations of metal ion and immobilized ligand, as well as GoBio Mini His-tag NTA Screening kits and GoBio Mini His-tag IDA Screening kits.

Target molecules

His-tagged proteins and other proteins with exposed histidine, cysteine and tryptophan on their surface.

See schematic depicting immobilized metal ion affinity chromatography.



Precharged IMAC resin resistant to DTT and EDTA




WorkBeads NiMAC

- Resin with extra strongly bound Ni²⁺ resulting in extremely low nickel ion leakage
- Highly resistant to reducing agents up to 20 mM DTT
- Highly resistant to chelating substances present in eukaryotic extracts or up to 20 mM EDTA
- Excellent purity, recovery and reproducible results
- Outstanding alkali stability with 0.5 M NaOH, extends the number of purification cycles
- Available in several different GoBio prepacked columns

Applications

Larger sample load incl. 20 mM DTT and 20 mM Na₂-EDTA

Column: GoBio Mini NiMAC 1 mL
 Sample: 50 mL His₆-GFP in binding buffer with 20 mM DTT and 20 mM Na₂-EDTA
 Binding buffer: 50 mM sodium phosphate, 300 mM NaCl, 10 mM imidazole, pH 8.0
 Elution buffer: 50 mM sodium phosphate, 300 mM NaCl, 300 mM imidazole, pH 8.0
 Elution: Step gradient, 100% elution buffer, 10 column volumes (CV)
 Flow rates: 0.5 mL/min (78 cm/h; elution); 1 mL/min (loading)



(A) Chromatogram with 50 mL load of His₆-GFP, (B) SDS-PAGE under reducing conditions of the feed, flow-through (FT) and eluted pool (1x: concentrated eluate, 5x: 1:5 diluted eluate) from 10 mL sample load and 50 mL sample load. (C) Comparison of target recovery and purity for the two different sample load purifications.

IMAC resins precharged with different metal ions



WorkBeads 40 Ni-NTA
WorkBeads 40 Co-NTA
WorkBeads 40 Zn-NTA
WorkBeads 40 Cu-NTA

WorkBeads 40 Ni-IDA
WorkBeads 40 Co-IDA
WorkBeads 40 Zn-IDA
WorkBeads 40 Cu-IDA

- Resins immobilized with either NTA (Nitrilotriacetic acid) or IDA (Iminodiacetic acid) and four different choices of metal ions Ni^{2+} , Co^{2+} , Zn^{2+} or Cu^{2+}
- Precharged with different metal ions for ease of use
- Low leakage of immobilized ligand and metal ions
- Resistant to harsh cleaning agents (NaOH). Note! The metal ions have to be stripped off before cleaning
- High binding capacity and flow rate
- Available in several different GoBio prepacked columns



Applications


Purification of His-tagged protein on GoBio Mini Ni-NTA vs. HisTrap™ FF +/- upstream GoBio Mini TREN

Purification of complex feeds expressed in different host cells can result in an extensive bioburden on the capture column in the form of DNA and different protein impurities when the feed is directly loaded without major pre-treatments. These impurities also often bind non-specifically to the target molecules and/or resin, and thus may be co-eluted with the final product.

One example of such a complex feed is His-tagged proteins expressed in bacteria. Therefore, WorkBeads 40 TREN prepacked in a GoBio Mini TREN column was introduced as a pre-treatment step upstream the IMAC column when purifying the *E. coli* translation initiation factor 3 (IF-3-His₆). Since this protein has a nucleic acid binding domain, host cell nucleic acids can potentially be a major co-eluting impurity in the eluates.

The GoBio Mini TREN column was operated in flow-through mode to capture impurities. This resin binds host cell nucleic acid (HCD), viruses and various host cell proteins (HCP), thereby reducing the foulant load on the subsequent IMAC column. The effect of the TREN column was evaluated by examining the removal of HCD and HCPs, see analysis.

Schematic view of GoBio Mini TREN upstream GoBio Mini Ni-NTA







IF-3-His₆ purification on GoBio Mini Ni-NTA and HisTrap FF +/- upstream GoBio Mini TREN

The effects of WorkBeads 40 TREN as a pre-treatment step:

- Increased purity in flow-throughs and eluates (Figures A-D)
- > 99.9% HCD removal in feed and 96-98% HCD removal in eluates
- No significant loss of target protein recovery (Figure B)
- ELISA showed that GoBio Mini TREN removes 49-62% more HCPs in eluates
- GoBio Mini Ni-NTA and HisTrap FF show analogous results (Figures B-D)

Analytical SEC of eluted IF-3-His₆ from only GoBio Mini Ni-NTA compared to GoBio Mini Ni-NTA in combination with GoBio Mini TREN

Sample volumes: 200 μ L
 Column: Superdex™ 200 10/300 GL (Cytiva)
 Buffer: 20 mM PBS, pH 7.4
 Flow rate: 30 cm/h (0.3 mL/min)



Conclusion

The final purity of the eluted His-tagged proteins, regarding HCD presence, was increased 97–98% when GoBio Mini TREN was placed upstream the IMAC column compared to a stand-alone step with only IMAC. Moreover, the purity was increased in terms of less co-eluted HCPs. In conclusion, the addition of an upstream pre-treatment step is highly advantageous to include early in a purification process to remove especially host-cell DNA that tend to interfere with downstream chromatographic processes.

Uncharged IMAC resins

WorkBeads 40 NTA WorkBeads 40 IDA


- Resins immobilized with either NTA (nitrilotriacetic acid) or IDA (iminodiacetic acid) for immobilization of your choice of metal ion
- Low leakage of immobilized ligand and metal ions of choice
- Resistant to harsh cleaning agents (NaOH). Note! The metal ions have to be stripped off before cleaning
- High binding capacity and flow rate
- Available in several different GoBio prepacked columns



Applications

Purification of clarified His₆-GFP on GoBio Mini Ni-NTA

Column: GoBio Mini Ni-NTA 1 mL
 Sample: 40 mL His₆-GFP in binding buffer
 Binding buffer: 50 mM sodium phosphate, 300 mM NaCl, 10 mM imidazole, pH 8.0
 Elution buffer: 50 mM sodium phosphate, 300 mM NaCl, 300 mM imidazole, pH 8.0
 Elution: 100% elution buffer in 5 CV
 Elution flow rate: 0.5 mL/min (75 cm/h)



(A) Chromatogram of the capture and elution of His₆-GFP. Absorbance at 280 nm (blue), absorbance at 490 nm (red) and percentage of elution buffer (green). (B) SDS-PAGE analysis of sample, flow-through, wash and eluted peak.



Technical specifications

WorkBeads NiMAC	
Target substance	His-tagged proteins
Matrix	Highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm
Precharged ions	Nickel (II)
Static binding capacity	> 80 mg/mL resin
Dynamic binding capacity ²	> 40 mg/mL resin
Metal ion capacity ³	> 60 $\mu\text{mol Cu}^{2+}$ /mL resin
Max flow rate (20 cm bed height and 5 bar)	600 cm/h
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, and additives such as 20 mM $\text{Na}_2\text{-EDTA}$, 20 mM dithiothreitol (DTT), 20 mM TECP, 20 mM β -mercaptoethanol, 8 M urea, 6 M guanidine-HCl, non-ionic detergents, 500 mM imidazole, 30% isopropanol, 0.5 M NaOH
pH stability	3 to 9 (working range) 2 to 14 (cleaning-in-place)
Storage	2 to 25°C in 20% ethanol

¹ Optimal flow rate during binding is depending on the sample. During column wash and elution, a flow rate of 1 mL/min and 5 mL/min can be used for 1 mL and 5 mL columns, respectively. Note: The maximum pressure the packed bed can withstand depends on the sample/liquid viscosity and chromatography resin characteristics. The pressure also depends on the tubing used to connect the column and the system restrictions after the column outlet.

² Maximum flow rate for aqueous buffers at 20°C. Decrease the maximum flow rate if the liquid has a higher viscosity. Higher viscosities can be caused by low temperature or presence of additives. Use half of the maximum flow rate for 20% ethanol for example.

	WorkBeads 40 Ni-NTA WorkBeads 40 Ni-IDA	WorkBeads 40 Co-NTA WorkBeads 40 Co-IDA	WorkBeads 40 Cu-NTA WorkBeads 40 Cu-IDA	WorkBeads 40 Zn-NTA WorkBeads 40 Zn-IDA
Matrix	Highly cross-linked agarose	Highly cross-linked agarose	Highly cross-linked agarose	Highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm	45 μm	45 μm
Chelating ligand	Nitrilotriacetic acid (NTA) or Iminodiacetic acid (IDA)	NTA or IDA	NTA or IDA	NTA or IDA
Metal ion	Nickel (II)	Cobalt (II)	Copper (II)	Zink (II)
Metal ion capacity for the chelating ligand ²	N/A	N/A	50 to 60 $\mu\text{mol Cu}^{2+}$ /mL (WorkBeads 40 Cu-IDA)	N/A
Dynamic binding capacity ³ (DBC)	> 60 mg His ₆ -GFP/mL resin	N/A	N/A	N/A
Maximum flow rate (20 cm bed height, 5 bar)	600 cm/h	600 cm/h	600 cm/h	600 cm/h
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, 8 M urea, 6 M guanidine-HCl, non-ionic detergents, 20% ethanol. Chelating substances (e.g. EDTA) will strip off the metal ions. Stripped resin: 10 mM HCl (pH 2), 10 mM NaOH (pH 12), 10 mM sodium citrate-HCl (pH 3).			
pH stability	7 to 9 (working) 2 to 12 (cleaning, stripped resin)	7 to 9 (working) 2 to 12 (cleaning, stripped resin)	7 to 9 (working) 2 to 12 (cleaning, stripped resin)	7 to 9 (working) 2 to 12 (cleaning, stripped resin)
Storage	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol	2 to 25°C in 20% ethanol

¹ The median particle size of the cumulative volume distribution.

² Metal ion capacity is determined by frontal analysis at 50% breakthrough using copper solution.

³ The binding capacity is determined using a GoBio Mini Ni-NTA 1 mL, equal value is expected for IDA resins. The binding capacity is dependent on the size of the target protein, and on the competition from impurities.

	WorkBeads 40 NTA	WorkBeads 40 IDA
Matrix	Highly cross-linked agarose	Highly cross-linked agarose
Average particle size ¹ (D_{v50})	45 μm	45 μm
Chelating ligand	Nitrilotriacetic acid (NTA)	Iminodiacetic acid (IDA)
Metal ion capacity ²	20 to 30 $\mu\text{mol Cu}^{2+}$ /mL resin	50 to 60 $\mu\text{mol Cu}^{2+}$ /mL resin
Maximum flow rate	600 cm/h (20 cm bed height, 5 bar)	600 cm/h (20 cm bed height, 5 bar)
Chemical stability	Compatible with all standard aqueous buffers used for protein purification, 8 M urea, 6 M guanidine-HCl, non-ionic detergents, 20% ethanol. Chelating substances (e.g. EDTA) will strip off the metal ions. Stripped resin: 10 mM HCl (pH 2), 10 mM NaOH (pH 12), 10 mM sodium citrate-HCl (pH 3).	
pH stability	7 to 9 (working range) 2 to 12 (cleaning, stripped resin)	7 to 9 (working range) 2 to 12 (cleaning, stripped resin)
Storage	2 to 25°C in ethanol	2 to 25°C in ethanol

¹ The median particle size of the cumulative volume distribution.

² Metal ion capacity is determined by frontal analysis at 50% breakthrough using copper solution.

Ordering information

Product name	Pack size	Article number
WorkBeads NiMAC	25 mL	40 653 001
	150 mL	40 653 003
	1 L	40 653 010
	5 L	40 653 050
WorkBeads 40 Ni-NTA	25 mL	40 651001
	150 mL	40 651003
	1 L	40 651010
	5 L	40 651050
	10 L	40 651060
WorkBeads 40 Co-NTA	25 mL	40 651401
	150 mL	40 651403
	1 L	40 651410
WorkBeads 40 Cu-NTA	25 mL	40 651301
	150 mL	40 651303
	1 L	40 651310
WorkBeads 40 Zn-NTA	25 mL	40 651501
	150 mL	40 651503
	1 L	40 651510
WorkBeads 40 Ni-IDA	25 mL	40 650 001
	150 mL	40 650 003
	1 L	40 650 010
	5 L	40 650 050
	10 L	40 650 060
WorkBeads 40 Co-IDA	25 mL	40 650 401
	150 mL	40 650 403
	1 L	40 650 410
WorkBeads 40 Cu-IDA	25 mL	40 650 301
	150 mL	40 650 303
	1 L	40 650 310
WorkBeads 40 Zn-IDA	25 mL	40 650 501
	150 mL	40 650 503
	1 L	40 650 510
WorkBeads 40 NTA	25 mL	40 602 001
	150 mL	40 602 003
	1 L	40 602 010
	5 L	40 602 050
	10 L	40 602 060
WorkBeads 40 IDA	25 mL	40 601001
	150 mL	40 601003
	1 L	40 601010
	5 L	40 601050
	10 L	40 601060

More information

Data Sheet, DS 40 653 010

WorkBeads NiMAC, GoBio prepacked columns

Data Sheet, DS 40 600 010

WorkBeads 40 NTA, WorkBeads 40 IDA, GoBio prepacked columns

Data Sheet, DS 40 650 010

WorkBeads 40 Ni-NTA, WorkBeads 40 Co-NTA, WorkBeads 40 Zn-NTA,
WorkBeads 40 Cu-NTA, GoBio prepacked columns

WorkBeads 40 Ni-IDA, WorkBeads 40 Co-IDA, WorkBeads 40 Zn-IDA,
WorkBeads 40 Cu-IDA, GoBio prepacked columns

→ [bio-works.com/product/imac-resin](https://www.bio-works.com/product/imac-resin)

