



2032 Concourse Drive San Jose, CA 95131
Re-order Phone: 888-813-7253
Re-order Fax: 402-452-5401
Send PO: orders@transgenomic.com
Tech Support: rjones@transgenomic.com

USER GUIDE

CARBOSep COREGEL 87C FAST ANALYSIS CARBOHYDRATE COLUMN

Ca⁺⁺ Form

CATALOG NUMBER CHO-99-5860

WARNING. THE TRANSGENOMIC CARBOSep COREGEL 87C FAST ANALYSIS CALCIUM FORM COLUMN IS PACKED WITH A POLYMERIC MATERIAL THAT REQUIRES SPECIAL CARE. INTRODUCTION OF ORGANIC SOLVENTS INTO THE COLUMN EXCEPT AS DESCRIBED BELOW WILL CAUSE THE POLYMER TO SWELL AND THE COLUMN WILL OVERPRESSURE. THE COLUMN WARRANTY IS INVALIDATED IF THIS OCCURS. CONSEQUENTLY, YOU SHOULD THOROUGHLY FLUSH THE LINES OF YOUR LC WITH ISO-PROPANOL FOLLOWED BY DE-IONIZED WATER TO REMOVE ANY ORGANIC SOLVENTS PRIOR TO INSTALLATION OF THIS COLUMN.

YOU SHOULD THOROUGHLY FAMILIARIZE YOURSELF WITH THE CONTENTS OF THIS MANUAL BEFORE USING YOUR COLUMN. IMPROPER USE WILL INVALIDATE THE WARRANTY.



2032 Concourse Drive San Jose, CA 95131
Re-order Phone: 888-813-7253
Re-order Fax: 402-452-5401
Send PO: orders@transgenomic.com
Tech Support: rjones@transgenomic.com

DESCRIPTION

The Transgenomic CARBOSep COREGEL 87C FAST ANALYSIS carbohydrate column contains a 7.8 x 100mm bed packed with a Transgenomic cation - exchange resin in the Ca⁺⁺ ionic form. It is specifically designed for the separation of polysaccharides, monosaccharides and sugar alcohols on a single chromatographic profile using water as the mobile phase. The primary mode of separation is steric exclusion although partitioning and ligand exchange mechanisms may be involved in the separation of certain monosaccharides. The column's physical dimensions (3/8 in. O.D. x 10 cm) are compatible with most commercially available heaters.

MOBILE PHASE

The only recommended mobile phase is water, de-gassed, de-ionized (especially free of metal ions and halides) and bacteria-free (filtered through 0.45 micron membrane). Alternatively, "HPLC grade" water is available commercially and is satisfactory for use with this column. The mobile phase should be de-gassed by vigorous heating prior to use and kept in a container that precludes introduction of airborne bacterial and fungal contamination. Fresh mobile phase should be prepared every 24 hours.

MOBILE PHASE FLOW RATE

It is good practice to limit mobile phase flow rates such that pump pressure does not exceed 68 bar (1000 psi). The recommended mobile phase flow rate for the CARBOSep COREGEL 87C FAST ANALYSIS column is 0.6 mL/min. A flow rate of 0.2 mL/min should be used during start-up until the column reaches the desired operating temperature. It is possible to employ mobile phase flow rates as high as 1.0 mL/min. Since all columns are tested at the factory prior to shipment, increased column pressure indicates that contaminants may have been deposited on the bed thereby restricting mobile phase flow.

MOBILE PHASE FLOW DIRECTION

An arrow is on the column label. This arrow is for reference purposes and indicates the preferred flow direction used during testing. However, the column can be operated with mobile phase flowing in either direction.

COLUMN PRESSURE

Remember that the pump pressure required to deliver mobile phases through the column is a *consequence* of mobile phase flow rate, column temperature, mobile phase viscosity,



2032 Concourse Drive San Jose, CA 95131
Re-order Phone: 888-813-7253
Re-order Fax: 402-452-5401
Send PO: orders@transgenomic.com
Tech Support: rjones@transgenomic.com

etc. Under normal operating conditions, a flow rate of 0.6 mL/min at 85° C should require pump pressures of less than 37 bar. It is inadvisable to utilize mobile phase flow rates that produce pump pressures in excess of 68 bar. If high pressures result from use of the column at nominal flow rates, this usually indicates that some contaminants have become deposited on the packing material and corrective action must be taken (see section on “Causes of Performance Loss”). To prevent irreversible damage to the column, however, you must exercise care in preparing mobile phases and samples. High column pressures nearly always results from improper use of the column. Use of a guard column (see below) will usually prevent contaminants from accumulating on the analytical column.

COLUMN TEMPERATURE

The Transgenomic CARBOSep COREGEL-87C FAST ANALYSIS carbohydrate column should always be used with a heating device. A particular characteristic of the packing material is the improved efficiency that results from use of the column at elevated temperatures. However, column temperature also influences sample retention and therefore must be carefully maintained to insure repeatable results. Although temperature can be used to influence certain separations, the temperature range of 85° - 90° C has been determined to be the optimum for the Transgenomic CARBOSep COREGEL 87C FAST ANALYSIS carbohydrate column. Temperatures below 80° C can be used for some applications but mobile phase pump pressures may be excessive unless flow rates are reduced. If it is necessary to use the column at lower than recommended temperatures, mobile phase flow rate should be adjusted to keep pump pressures below 68 bar (1000 psi).

PRE-COLUMN FILTER

Pre-column filters containing 0.5 – 2.0 micron porosity passivated stainless steel or titanium frits should be used between the sample injector and the column to remove particulates from the mobile phase stream. This will help prevent excessive pressure from developing through the analytical column and will prolong column life.

GUARD SYSTEMS

Guard systems should be used with polymeric carbohydrate columns because sample and mobile phase contamination can result in excessive column pressures. Contaminants such as salts and proteins can alter column performance and should always be removed from samples prior to injection onto the column. We recommend the Transgenomic CARBOSep COREGEL-87C Guard Kit (PN: CHO-99-2360). This is a guard cartridge system that contains the same polymer used in the CARBOSep COREGEL-87C



2032 Concourse Drive San Jose, CA 95131
Re-order Phone: 888-813-7253
Re-order Fax: 402-452-5401
Send PO: orders@transgenomic.com
Tech Support: rjones@transgenomic.com

analytical column. The Guard Kit includes a stainless steel Guard Cartridge Holder and a 2/pk of 87C Guard Cartridges (PEEK, 4.0x20mm). Replacement 87C Guard Cartridges (2/pk) are available (PN: CHO-99-1360). Cartridge replacement is required when increased column pressure and/or loss of resolution is observed. Silica guard columns are not recommended due to degradation and eventual leakage into the analytical column. An alternative with more capacity is the Transgenomic CARBOSep COREGEL-87C Guard Column (PN: CHO-99-3560) which is 4.6 x 50mm. Use of guard columns should dramatically extend column lifetime and column cleaning or regeneration (described below) should not be required.

SAMPLE PREPARATION

The key to long column life is proper treatment of sample prior to injection onto the polymer bed. You should avoid introduction of fats, oils, proteinaceous materials and heavy metal ions into the column by either mobile phases or samples. If possible, you should avoid introduction of particulate matter onto the column. These will ultimately cause an increase in operating pressure and may be difficult or impossible to remove. We recommend sample treatment in order to prolong column life.

SAMPLE VOLUME

The Transgenomic CARBOSep COREGEL-87C FAST ANALYSIS carbohydrate column contains a bed of 0.78 x 10-cm. Although this bed should tolerate large injection volumes, remember that one of the separation modes is steric exclusion. Under these circumstances, smaller sample volumes usually promote higher separation efficiencies. Consequently, you must determine empirically the maximum injection volume tolerated by the column for your particular sample. Although we generally use sample volumes in the 10 – 50 μ L range, problems should not occur with sample volumes of as high as 100 μ L. Injections of 500 μ L or more may cause problems, depending on sample content. If the injection volume is too great, peaks may broaden or merge with nearby peaks.

DETECTION AND SENSITIVITY

The mobile phase requirement for Transgenomic CARBOSep COREGEL-87C carbohydrate column enables chromatographers to use a variety of detectors. Spectrophotometers, refractometers and electrochemical detectors can all be used successfully. If electrochemical detectors are used, note that high temperatures may be incompatible with some working electrodes. Selection of the right electrode with desired



2032 Concourse Drive San Jose, CA 95131
Re-order Phone: 888-813-7253
Re-order Fax: 402-452-5401
Send PO: orders@transgenomic.com
Tech Support: rjones@transgenomic.com

temperature and mobile phase strength is the key to successful applications. If higher sensitivity is required, post-column reactions followed by the appropriate detector for the reaction product (e.g. fluorometer, photometer) can be utilized. Remember that sensitivity of detection is ultimately determined by the type of detector chosen; the responsibility of the column is simply to separate the compounds of interest.

STORAGE

The column as supplied is equilibrated with de-ionized water. This is also the recommended mobile phase for storage. When the column is stored, be sure the end fittings are tightly sealed using the end plugs provided. Long term storage without these seals can result in partial drying of the packing material and high pressures can ensue. Under these circumstances invert the column and pump de-ionized water at a flow rate of 0.1 mL/min at 90° C. Gradually increase the flow rate to 0.5 mL/min. Normal pressure should be observed and the column can be used in either direction. If this does not correct the problem, the column may have become contaminated with particulates or other material.

POSSIBLE CAUSES OF PERFORMANCE LOSS IN THE TRANSGENOMIC CARBOSep COREGEL-87C CARBOHYDRATE FAST ANALYSIS CALCIUM FORM COLUMN

The following outline is intended as an aid in locating sources of performance loss. Because of the nature of polymeric materials and the manufacturing procedures employed by Transgenomic, it is highly unusual for a column to lose performance due to manufacturing problems. In our experience, nearly all column failures are a result of introduction of contaminants onto the polymer bed. Use of a guard column will help prevent these problems, as will sample pretreatment (see above). All Transgenomic columns are thoroughly tested prior to shipment and are supplied with a sample chromatogram illustrating performance of the particular column.

1. Post-column mixing and or diffusion – make sure the tubing length and tubing I.D. are kept to a minimum (ref. 2).
2. Improper column temperature
3. Improper mobile phase flow rate



2032 Concourse Drive San Jose, CA 95131
Re-order Phone: 888-813-7253
Re-order Fax: 402-452-5401
Send PO: orders@transgenomic.com
Tech Support: rjones@transgenomic.com

4. Insufficient equilibration time with mobile phase
5. Improper pH of mobile phase
6. Presence of cations in eluent – (e.g., Na⁺ or H⁺)
7. Polymer contamination
 - a. High column pressure accompanies performance loss
 - (1) particulate accumulation on inlet frit or polymer bed
 - (a) sample origin — filter or centrifuge
 - (b) mobile phase origin — filter mobile phase; enclose mobile phase reservoirs
 - (c) system origin — flush all lines and pump; install in-line system filter
 - (2) proteinaceous material accumulation
 - (a) microbial growth in samples
 - (b) microbial growth in mobile phase
 - (3) refractory organic component from sweetener samples
 - b. Normal column pressure accompanies performance loss
 - (1) metal ion contamination
 - (a) inappropriate steel alloy present in LC system
 - (b) samples contaminated with metal ions
 - (2) organic contamination
 - (a) fats, oils, lipids in sample – polymer surface becomes coated
 - (b) non-specific organic from improperly prepared mobile phase or source material
8. Bed compression (voids)
 - a. Excessive mobile phase flow rate
 - b. Use of organic modifier (not recommended)

**OPERATIONS DESIGNED TO CORRECT PERFORMANCE LOSSES
RESULTING FROM POLYMER CONTAMINATION OR BED COMPRESSION**

1. **Prepare Fresh Mobile Phase.** In some cases, performance loss is traced to mobile phase contamination. Therefore, prepare fresh mobile phase and flush all liquid lines before using column; mobile phase should be filtered through 0.2 to 0.4 micron membranes and de-gassed prior to use.
2. **“Loosening” the Polymer bed.** Many polymers lack the rigidity associated with silica materials and can compress or collapse if inappropriately high mobile phase flow rates are utilized. They are resilient, however, and the compression is reversible except in severe cases. To correct collapsed beds, shut off the pump and allow the polymer to “relax” for approximately 30 minutes. Invert the column and mobile phase at 0.1mL/min overnight at 90° C. Return column to normal operating conditions.
3. **Cleaning of Polymer.** If performance problems persist and particularly if high column pressures remain, an attempt should be made to clean the column to remove built up contaminants. Cleaning and regeneration procedures outlined below will in some cases restore performance by removing contaminants from the polymer bed. It is important, however, to attempt to locate the source of the problem before again using the column for analysis of samples.

Prepare an aqueous solution of 20% acetonitrile. Set column temperature to 65° C and pump solution through (inverted) column at 0.1 mL/min overnight. If necessary, adjust mobile phase flow rate such that pressure does not exceed 55 bar. You may see a dark – colored material eluting from the column.

Next day, replace acetonitrile solution with de-ionized water, and continue pumping at 0.1 mL/min to determine if high pressure has subsided. If pressure is low, return column temperature to 90 °C and gradually increase mobile phase flow rate to 0.6 ml/min. Test column under normal but performance remains inadequate, attempt regeneration procedure described below. If pressure does not return to normal, column may be permanently damaged and require replacement.

4. **Column Regeneration.** Prepare an aqueous 1% solution of calcium nitrate containing 0.25% EDTA. Filter this solution through a 0.45 micron membrane and degas. With the column in the inverted direction, pump at 0.1 ml/min at 90° C overnight. Retest the column under normal operating conditions. Note that it may take some time for the baseline to stabilize.
5. **Column Replacement.** Above procedures will restore performance only in certain cases. Heavy metal contamination and certain organic contaminants are particularly refractory and may not respond to treatment. Under these



2032 Concourse Drive San Jose, CA 95131
Re-order Phone: 888-813-7253
Re-order Fax: 402-452-5401
Send PO: orders@transgenomic.com
Tech Support: rjones@transgenomic.com

circumstances, column replacement is necessary. It is highly advisable to locate the cause of the problem before installing a new column. Consult the manufacturer of your LC system for aid in this matter.

WARRANTY

Transgenomic Incorporation warrants its columns to perform within specifications for 30 days. Warranty does not cover acts of neglect such as, but not limited to: use of improper or improperly prepared eluents; overpressure resulting from excessive eluent flow rate; overpressure or loss in performance resulting from improperly prepared samples; loss in performance resulting from improper column storage or cleaning procedure.

REFERENCES

- (1) Benson, J.R. and Woo, D.J. (1984) *J. Chrom. Sci.* 22 (9) 386-399
- (2) Majors, R. (1983), *LC* 1 (Nov), 464

2013 Transgenomic, Inc. All rights reserved. Reproduction in whole or in part is strictly prohibited.
