







# 1. Reversed Phase Specialty Columns


## Specifications

Packing Material	Cholester			PBr			PFP	$\pi$ NAP			PYE	NPE		
	Fully-Porous	Core-Shell	Core-Shell	Fully-Porous	Core-Shell	Core-Shell	Fully-Porous	Fully-Porous	Fully-Porous	Fully-Porous	Fully-Porous	Fully-Porous		
Average Particle Size ( $\mu\text{m}$ )	2.5	3	5	2.6	3	5	2.6	5	2.5	3	5	5	5	
Average Pore Size ( $\text{\AA}$ )	130	120	90	120	120	90	120	120	130	120	120	120	120	
Specific Surface Area ( $\text{m}^2/\text{g}$ )	330	300	150	300	300	150	300	300	330	300	300	300	300	
Bonded Phase Structure														
Bonded Phase	Cholesteryl group			Pentabromobenzyl group			Pentafluorophenyl group	Naphthylethyl group	Pyrenylethyl group	Nitrophenylethyl group				
Bonding Type	Monomeric													
Main Interaction	Hydrophobic interaction Molecular shape selectivity			Hydrophobic interaction Dispersion force			Hydrophobic interaction $\pi$ - $\pi$ interaction Dipole-dipole interaction	Hydrophobic interaction $\pi$ - $\pi$ interaction	Hydrophobic interaction $\pi$ - $\pi$ interaction Dispersion force Molecular shape selectivity	Hydrophobic interaction $\pi$ - $\pi$ interaction Dipole-dipole interaction				
End-Capping	Near-perfect treatment													
Carbon Content	21%	20%	—	8%	—	—	10%	14%	11%	18%	18%	9%	9%	
Usable pH Range	2 ~ 7.5													
Features	<ul style="list-style-type: none"> <li>Usable under the same conditions as <math>\text{C}_{18}</math></li> <li>High molecular shape selectivity</li> </ul>			<ul style="list-style-type: none"> <li>Separate hydrophilic compounds under reversed-phase conditions</li> <li>Separate using dispersion force</li> </ul>			<ul style="list-style-type: none"> <li>Weak dipole-dipole interaction</li> </ul>	<ul style="list-style-type: none"> <li>Stronger <math>\pi</math>-<math>\pi</math> interaction than phenyl column</li> </ul>	<ul style="list-style-type: none"> <li>Very strong <math>\pi</math>-<math>\pi</math> interaction</li> </ul>	<ul style="list-style-type: none"> <li>Strong dipole-dipole interaction</li> </ul>				

\* Silica Gel : Fully-Porous...High purity porous spherical silica    Core-Shell...Core-Shell silica gel

## Selectivity for positional isomers of dinitrobenzene

Different stationary phases exhibit different selectivity due to the use of forces that  $\text{C}_{18}$  (hydrophobic interaction) does not have. By using these columns, you can achieve separation that cannot be done using only  $\text{C}_{18}$ .



# COSMOSIL Cholester / COSMOCORE Cholester




- Cholesterol-bonded stationary phase
- Increased stereoselectivity and improved resolution for geometric isomers
- Usable under the same conditions as C<sub>18</sub>

## Suitable Samples

- Natural compounds
- Structurally similar compounds
- Polyphenols, catechins, fat-soluble vitamins and flavones

## Hydrophobic Interaction

The below figure shows the comparison of hydrophobic interactions with competitor C<sub>18</sub> columns. Cholester provides about the same hydrophobicity as alkyl group-bonded types (C<sub>18</sub>, C<sub>30</sub>). It is not necessary to change the analytical conditions when replacing C<sub>18</sub> or C<sub>30</sub> columns with Cholester.



## Molecular Shape Selectivity

The stationary phase of Cholester has a very rigid structure and can distinguish different molecular shapes. Cholester retains planar triphenylene longer than non-planar *o*-Terphenyl.






Figure. Comparison of molecular shape selectivity



## Improvement in Separation


COSMOSIL Cholester provides enhanced selectivity over traditional C<sub>18</sub> columns and offers greater performance in separating isomers or other closely related compounds. COSMOSIL Cholester is ideal for method development and serves as an excellent alternative to traditional C<sub>18</sub> columns. The figure below shows analytical data of chloroacetophenone isomers. These isomers are difficult to separate with C<sub>18</sub> and C<sub>30</sub>, but they are well resolved by COSMOSIL Cholester.

### ● Comparison with competitor's C<sub>18</sub> and C<sub>30</sub> columns



## Efficiency of Preparative Separation


The figure below shows the comparison of efficiency of preparative separation with a C<sub>18</sub> column. Both columns show good separation. However, sample loading capacity for preparative separations can be affected by a slight difference in separation ability. COSMOSIL Cholester can load 4 times the sample volume compared with C<sub>18</sub> columns.




## About Core-Shell 2.6Cholester Particles

COSMOCORE 2.6Cholester is packed with cholesterol-bonded 2.6 µm core-shell particles. It delivers performance equivalent to sub-2 µm particles at faster flow rate and analysis time while maintaining a lower back pressure. COSMOCORE can also be used in longer column size to gain additional resolution.

### ● Catechins (Standard)




### ● Comparison of Column Pressure




Column Size : 2.1 mmI.D. x 100 mm  
 Mobile Phase : Acetonitrile/Water = 70/30  
 Temperature : 40°C

## Applications


### ● Catechins (Commercial Green Tea)




### ● Carotene



### ● Fat-Soluble Vitamins




### ● Flavone




## Applications

### ● Saikosaponins



### ● Tomato Extracts



## Ordering Information

### ● COSMOSIL Cholester Analytical / Preparative Columns (Particle Size: 5 µm)

#### Packed Column

I.D. x Length (mm)	Product Number
1.0 × 150	05968-71
1.0 × 250	05969-61
2.0 × 30	08565-51
2.0 × 50	06352-91
2.0 × 100	06948-01
2.0 × 150	05971-11
2.0 × 250	05972-01
3.0 × 150	05973-91
3.0 × 250	05974-81
4.6 × 50	06359-21
4.6 × 100	06591-61

I.D. x Length (mm)	Product Number
4.6 × 150*	05976-61
4.6 × 150 3 lots set*	07970-03
4.6 × 250*	05977-51
10 × 50	16590-61
10 × 100	16591-51
10 × 150	08011-91
10 × 250	05979-31
20 × 50	05981-81
20 × 100	15995-01
20 × 150	06088-71
20 × 250	05982-71

I.D. x Length (mm)	Product Number
28 × 100	16592-41
28 × 150	16593-31
28 × 250	05985-41

#### Guard Column / Guard Cartridge

I.D. x Length (mm)	Product Number
4.6 × 10	05975-71
4.6 × 10 Cartridge †	19183-24
10 × 20	05978-41
20 × 20	05980-91
28 × 50	05983-61

\* Columns for validation

† 2 cartridges included. Guard cartridge holder required; refer to page 78

### ● COSMOSIL 3Cholester Fast LC Columns (Particle Size: 3 µm)

#### Packed Column

I.D. x Length (mm)	Product Number
2.0 × 50	19188-61
2.0 × 75	19189-51
2.0 × 100	19190-11
2.0 × 150	19191-01
2.0 × 250	19192-91
3.0 × 50	19194-71
3.0 × 75	19195-61

I.D. x Length (mm)	Product Number
3.0 × 100	19196-51
3.0 × 150	19197-41
3.0 × 250	19198-31
4.6 × 50	19199-21
4.6 × 75	19200-71
4.6 × 100	19300-61
4.6 × 150	19151-21
4.6 × 250	19316-71

#### Guard Cartridge

I.D. x Length (mm)	Product Number
2.0 × 10 Cartridge*	19325-64
4.6 × 10 Cartridge*	19344-14

\* 2 cartridges included. Guard cartridge holder required; refer to page 78.

### ● COSMOSIL 2.5Cholester Fast LC columns (Particle Size: 2.5 µm)

#### Packed Column

I.D. x Length (mm)	Product Number
2.0 × 50	09000-01
2.0 × 75	09047-11

I.D. x Length (mm)	Product Number
2.0 × 100	09048-01
3.0 × 50	09049-91

I.D. x Length (mm)	Product Number
3.0 × 75	09050-51
3.0 × 100	09051-41

### ● COSMOCORE 2.6Cholester UHPLC Columns (Particle Size: 2.6 µm)

#### Packed Column

I.D. x Length (mm)	Product Number
2.1 × 30	12858-91
2.1 × 50	12859-81
2.1 × 75	12860-41
2.1 × 100	12861-31
2.1 × 150	12862-21

I.D. x Length (mm)	Product Number
3.0 × 30	12863-11
3.0 × 50	12864-01
3.0 × 75	12866-81
3.0 × 100	12867-71
3.0 × 150	12868-61

I.D. x Length (mm)	Product Number
4.6 × 30	12869-51
4.6 × 50	12870-11
4.6 × 75	12871-01
4.6 × 100	12872-91
4.6 × 150	12873-81
4.6 × 250	12875-61

COSMOCORE's connector is the same type as Waters UPLC® columns.