

ACE® 300Å HPLC Columns

Bonded Phase	USP Listing	Particle Size (µm)	Pore Size (Å)	pH Range	Temp. Max °C (low pH)	Temp. Max °C (high pH)	Carbon Load (%)	SurfaceArea (m2/g)	Endcapped
C18-300	L1	3 5	300	2-8	60	40	9	100	YES
C8-300	L7	3 5	300	2-8	60	40	5	100	YES
C4-300	L26	3 5	300	2-8	60	40	2.6	100	YES
Phenyl-300	L11	3 5	300	2-8	60	40	5.3	100	YES
CN-300	L10	3 5	300	2-7	60	40	2.6	100	YES