

## UltraCore 產品規格

| Bonded Phase      | USP Listing | Particle Size (µm) | Pore Size (Å) | pH Range | Temp. Max °C (low pH) | Temp. Max °C (high pH) | Carbon Load (%) | Surface Area (m2/g) | Endcapped    |
|-------------------|-------------|--------------------|---------------|----------|-----------------------|------------------------|-----------------|---------------------|--------------|
| SuperC18          | L1          | 2.5<br>5           | 95            | 1.5-11.0 | 60                    | 40                     | 7<br>5.4        | 130<br>100          | Encapsulated |
| Super PhenylHexyl | L11         | 2.5<br>5           | 95            | 1.5-11.0 | 60                    | 40                     | 4.6<br>3.6      | 130<br>100          | Encapsulated |
| C18               | L1          | 3.5                | 95            | 2-9      | 60                    | 40                     | 8.2             | 115                 | YES          |
| Phenylhexyl       | L11         | 3.5                | 95            | 2-9      | 60                    | 40                     | 6.2             | 115                 | YES          |
| Biphenyl          | L11         | 3.5                | 95            | 2-9      | 60                    | 40                     | 6.5             | 115                 | YES          |
| C18-Amide         | L60         | 3.5                | 95            | 2-9      | 60                    | 40                     | 5.5             | 120                 | YES          |