

ZORBAX 300 StableBond (300SB) Agilent 產品規格表

Bonded Phase	Pore Size	Surface Area (m ² /g)	Temp. Limits	pH Range	End capped	CarbonLoad	USP
ZORBAX 300SB-C18	300	45	90 °C	1.0 - 8.0	No	2.80%	L1
ZORBAX 300SB-C8	300	45	80 °C	1.0 - 8.0	No	1.50%	L7
ZORBAX 300SB-C3	300	45	80 °C	1.0 - 8.0	No	1.10%	L56
ZORBAX 300SB-CN	300	45	80 °C	1.0 - 8.0	No	1.20%	L10
ZORBAX 300-Diphenyl	300	45	80 °C	1.0 - 8.0	Yes	1.90%	L11