

ZORBAX Eclipse Plus Agilent 產品規格表

Bonded Phase	Pore Size	Surface Area (m ² /g)	Temp. Limits	pH Range	End capped	Carbon Load	USP
ZORBAX Eclipse Plus C18	95Å	160	60 °C	2.0 - 9.0	Double	9%	L1
ZORBAX Eclipse Plus C8	95Å	160	60 °C	2.0 - 9.0	Double	7%	L7
ZORBAX Eclipse Plus Phenyl-Hexyl	95Å	160	60 °C	2.0 - 8.0	Double	9%	L11
ZORBAX Eclipse PAH	95Å	160	60 °C	2.0 - 8.0	No	14%	L1