

ZORBAX Eclipse XDB Agilent 產品規格表

Bonded Phase	Pore Size	Surface Area (m ² /g)	Temp. Limits	pH Range	End capped	Carbon Load	USP
ZORBAX Eclipse XDB C18	80	180	60 °C	2.0 - 9.0	Double	10%	L1
ZORBAX Eclipse XDB C8	80	180	60 °C	2.0 - 9.0	Double	7.60%	L7
ZORBAX Eclipse XDB Phenyl	80	180	60 °C	2.0 - 9.0	Double	7.20%	L11
ZORBAX XDB CN	80	180	60 °C	2.0 - 8.0	Double	4.30%	L10