

ZORBAX GF-250 and GF-450 Gel Filtration columns Agilent 產品規格表

Bonded Phase	Pore Size	Particle Size (μm)	MW Range	Surface Area (m ² /g)	pH Range	FlowRate	Max Pressure	USP
ZORBAX GF-250	150	4	4,000-400,000	140	3.0 - 8.0	<3.0 mL/min	350 bar	L33 L35
ZORBAX GF-450	300	6	10,000-900,000	50	3.0 - 8.0	<3.0 mL/min	350 bar	L35