

ZORBAX StableBond (SB) Agilent 產品規格表

Bonded Phase	Pore Size	Surface Area (m ² /g)	Temp. Limits	pH Range	End capped	CarbonLoad	USP
ZORBAX SB C18	80	180	60 °C	0.8 - 8.0	No	10%	L1
ZORBAX SB C8	80	180	60 °C	1.0 - 8.0	No	5.50%	L7
ZORBAX SB C3	80	180	60 °C	1.0 - 8.0	No	4%	L56
ZORBAX SB Phenyl	80	180	60 °C	1.0 - 8.0	No	5.50%	L11
ZORBAX SB CN	80	180	60 °C	1.0 - 8.0	No	4%	L10
ZORBAX SB CN	80	180	60 °C	1.0 - 8.0	No	proprietary	