

Bioseparation MCI 三菱
Ion exchange chromatography: MCI GEL™ CQK and CQA series
產品規格表

Product Name	Base Material	Functional group	Counter ion	Mean particle size (μm)	Pore size (nm)	pH range	Typical Application
CQA31P	HMA	DEAE	Cl ⁻	30	60	<11	Proteins
CQA35P	HMA	QA	Cl ⁻	30	60	<11	Proteins
CQK30P	HMA	SP	NA ⁺	30	60	<11	Proteins
CQK31P	HMA	CM	NA ⁺	30	60	<11	Proteins