















Reverse Phase 逆相型 Daicel 產品規格表

 Amylose derivative coated on silica gel R: 	CHIRALPAK® AD-3R / CHIRALPAK® AD-RH USP L51 AD為多糖衍生物逆相塗敷型管柱，矽膠表面塗敷有直鏈澱粉-三（3,5-二甲基苯基氨基甲酸酯）
 Amylose derivative coated on silica gel R: 	CHIRALPAK® AS-3R / CHIRALPAK® AS-RH AS為多糖衍生物逆相塗敷型管柱，矽膠表面塗敷有直鏈澱粉-三（S）-α-甲基苯基氨基甲酸酯）
 Amylose derivative coated on silica gel R: 	CHIRALPAK® AY-3R / CHIRALPAK® AY-RH AY為多糖衍生物逆相塗敷型管柱，矽膠表面塗敷有直鏈澱粉-三（5-氯-2-甲基苯基氨基甲酸酯）
 Amylose derivative coated on silica gel R: 	CHIRALPAK® AZ-3R / CHIRALPAK® AZ-RH AZ為多糖衍生物逆相塗敷型管柱，矽膠表面塗敷有直鏈澱粉（3-氯-4-甲基苯基氨基甲酸酯）
 Cellulose derivative coated on silica gel R: 	CHIRALCEL® OD-3R / CHIRALCEL® OD-RH / CHIRALCEL® OD-R USP L4 OD為多糖衍生物逆相塗敷型管柱，矽膠表面塗敷有纖維素-三（3,5-二甲基苯基氨基甲酸酯）
 Cellulose derivative coated on silica gel R: 	CHIRALCEL® OJ-3R / CHIRALCEL® OJ-RH USP L80 OJ為多糖衍生物逆相塗敷型管柱，矽膠表面塗敷有纖維素-三（4-甲基苯甲酸酯）
 Cellulose derivative coated on silica gel R: 	CHIRALCEL® OZ-3R / CHIRALCEL® OZ-RH OZ為多糖衍生物逆相塗敷型管柱，矽膠表面塗敷有纖維素-三（3-氯-4-甲基苯基氨基甲酸酯）
 Cellulose derivative coated on silica gel R: 	CHIRALCEL® OX-3R / CHIRALCEL® OX-RH OX為多糖衍生物逆相塗敷型管柱，矽膠表面塗敷有纖維素-三（4-氯-3-甲基苯基氨基甲酸酯）