








Special Applications 特殊應用 Daicel 產品規格表

	CHIRALPAK® ZWIX(+) / CHIRALPAK® ZWIX(-) ZWIX(+)矽膠表面共價鍵合-奎寧(8S,9R)-(1S,2S)-環己基氨基磺酸衍生物 ZWIX(-)矽膠表面共價鍵合-奎寧(8R,9S)-(1R,2R)-環己基氨基磺酸衍生物
CROWNPAK® CR-I(+)  CROWNPAK® CR-I(-) 	CROWNPAK® CR-I(+) / CROWNPAK® CR-I(-) CR-I矽膠表面塗覆鏡像冠醚
 CR (+) CR (-)	ACROWNPAK® CR(+) / CROWNPAK® CR(-) CR(+) USP L66 CR矽膠表面塗覆鏡像冠醚
	CHIRALPAK® QD-AX / CHIRALPAK® QN-AX QD-AX 矽膠表面共價鍵合-O-9-(叔丁酯氨基甲酯)奎寧定
	CHIRALPAK® WH USP L32 WH矽膠表面共價鍵合有L-脯氨酸衍生物
	CHIRALPAK® MA(+) MA矽膠表面塗敷有L-丙氨酸衍生物
	CHIRALPAK® OP(+) OP矽膠表面塗敷-二苯基-(2-吡啶)甲基-甲基丙烯酸酯
	CHIRALPAK® OT(+) OT矽膠表面塗敷-三苯基-甲基-甲基丙烯酸酯