

Davisil Media Grace 產品規格表

Normal Pore Size									
Characteristic	35Å	60Å	150Å	250Å	500Å	1000Å	1500Å	2500Å	4500Å
Surface Area (m ² /g)	700	550	330	285	80	40	25	17	10
Pore Volume (mL / g)	0.6	0.9	1.2	1.8	1.1	1.1	1.1	1.1	1.1
pH (5% suspension)	5.7	7.3	7.3	7.5	8	9	9	9	9.5
H ₂ O (weight %) [†]	< 6%	< 6%	< 6%	< 6%	< 6%	< 6%	< 6%	< 6%	< 6%
Bulk Density (kg / m ³)	720	530	350	210	370	370	370	370	370