

Allsphere 系列 Dr. Maisch 品規格表

Phase	Particle Size (μm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Loading	Endcapping	USP
ODS-1	5	80	220	7%	Partial	L1
ODS-2	3 5	80	220	12%	YES	L1
C8	3 5	80	220	6%	YES	L7
C6	5	80	220	4%	YES	L15
C1 (TMS)	5	80	220	3%	NO	L13
Phenyl	5	80	220	3%	YES	L11
Cyano	5	80	220	3.50%	NO	L10
Amino (NH ₂)	5	80	220	3%	NO	L8
Silica	3 5	80	220	-	NO	L3
SAX	5	100	220	4%	NO	-
SCX	5	100	220	4%	NO	-