

Brava 系列 Dr. Maisch 產品規格表

Phase	Particle Size (μm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Loading	Endcapping	USP
C18 BDS	3 5	145	185	8.50%	YES	L1
C18 ODS	3 5	130	195	8.50%	YES	L1
C8	3 5	130	195	6%	YES	L7
C8 BDS	3 5	145	185	5.50%	YES	L7
Phenyl	5	130	195	-	NO	L11
Cyano	5	130	195	-	NO	L10
Cyano BDS	5	145	185	-	NO	L10
Amino (NH ₂)	5	130	195	-	NO	L8
Silica	5	130	195	-	NO	L3