

Grom Sil 系列 Dr. Maisch 產品規格表

Phase	Particle Size (µm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Loading	Endcapping	USP
ODS-3 CP	3	120 300	320 170	15% 6%	NO	L1
	5					
	7					
	10					
ODS-4 HE	3	120 200	300 200	16% 11%	YES	L1
	4					
	5					
	7 10					
ODS-5 ST	3	60 120 200	580 300 200	22% 17% 12%	YES	L1
	4					
	5					
	7 10					
ODS-7 pH	4	80	510	22%	NO	L1
Octyl-5 CP	3	120 300	320 170	10% 5.5%	NO	L137
	5					
	7					
	10					
Phenyl-1 FE	3	120 300	320 150	9.50%	YES	L11
	5					
	10					
Phenyl-2 CP	5	120 300	320 170	7.40%	NO	L11
Phenyl-3 CP	5	120	320	4%	NO	-
Amino-2 PA	5	120	300	-	NO	L8
Diol	5	60	580	-	NO	L20
	10	120 200	300 200			