

Platinum 系列 Dr. Maisch 產品規格

Phase	Particle Size (μm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Loading	Endcapping	USP
C18	1.5 3 5	100	200	6%	YES	L1
EPS C18	1.5 3 5	100	200	5%	NO	L1
C8	1.5 3 5	100	200	4%	YES	L7
EPS C8	3 5	100	200	2.50%	NO	L7
Phenyl	3 5	100	200	-	YES	L11
Cyano	3 5	100	200	-	NO	L10
Amino (NH ₂)	3 5	100	200	-	NO	L8
Silica	3 5	100	200	-	-	L3
SAX	3 5	100	200	-	NO	-