

VisionHT 系列 Dr. Maisch 產品規格表

Phase	Particle Size (μm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Loading	Endcapping	USP
C18 HighLoad	1.5 3 5 10	120	220	11%	YES	L1
C18 Basic	1.5 3 5 10	120	220	5%	Proprietary	L1
C18 Classic	1.5 3 5 10	100	200	6%	YES	L1
C18 Polar	1.5 3 5 10	100	200	5%	NO	L1
HILIC	1.5 3 5 10	120	220	-	NO	L3
Silica	1.5 3 5 10	120	220	-	NO	L3