

Alltima 系列 Avantor Hichrom 產品規格表

Phase	Particle Size (μm)	Pore Size ($\text{\AA}$)	Surface Area (m^2/g)	Carbon Loading	Endcapping	USP
C18	3 5 10	100	340	16%	YES	L1
C18 LL	5	100	340	9%	YES	L1
C8	3 5 10	100	340	9%	YES	L7
Phenyl	3 5	100	340	7.5%	YES	L11
Cyano	3 5	100	340	-	YES	L10
Amino	3 5	100	340	-	NO	L8
Silica	3 5 10	100	340	-	NO	L3