

Alltima HP 系列 Avantor Hichrom 產品規格表

Phase	Particle Size (μm)	Pore Size ($\text{\AA}$)	Surface Area (m^2/g)	Carbon Loading	Endcapping	USP
C18	3 5	190	200	12%	YES	L1
C18 EPS	3 5	190	200	4%	YES	L1
C18 Hi-Load	3 5	100	450	24%	YES	L1
C18 AQ	3 5	100	450	20%	YES	L1
C18 Amide	3 5	190	200	12%	YES	L1
C8	3 5	190	200	8%	YES	L7
Cyano	3 5	190	200	4%	YES	L10
Silica	3 5	100	450	-	NO	L3
HILIC	1.5 3 5	120	230	-	NO	L3