

Genesis 系列 Avantor Hichrom 產品規格表

Phase	Particle Size (μm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Loading	Endcapping	USP
C18	3 4 7	120	300	18%	YES	L1
C8	3 4	120	300	11%	NO	L7
C8 e/c	4	120	300	11%	YES	L7
AQ	4	120	300	15%	YES	L1
Phenyl	4	120	300	9.4%	YES	L11
CN	4	120	300	7%	YES	L10
NH2	4	120	300	3.5%	NO	L8
Silica	4	120	300	-	-	L3