

Prevail 系列 Avantor Hichrom 產品規格表

Phase	Particle Size (μm)	Pore Size ($\text{\AA}$)	Surface Area (m^2/g)	Carbon Loading	Endcapping	USP
C18 Select	3 5	110	350	17%	YES	L1
C18	3 5	110	350	15%	YES	L1
C8	3 5	110	350	8%	YES	L7
Phenyl	3 5	110	350	7%	YES	L11
Cyano	3 5	110	350	-	YES	L10
Amino	3 5	110	350	-	NO	L8
Silica	3 5	110	350	-	-	L3
Organic Acid	3 5	110	350	-	YES	-