

Vydac MS 系列 Avantor Hichrom 產品規格表

Phase	Particle Size (μm)	Pore Size ($\text{\AA}$)	Surface Area (m^2/g)	Carbon Loading	Endcapping	USP
218MS C18	3 5 10	300	60-110	8%	YES	L1
238MS C18	5	300	70	4%	YES	L1
208MS C8	5	300	70	5%	YES	L7
214MS C4	5 10	300	70-110	3%	YES	L26
219MS Diphenyl	5	300	70	4%	YES	L11