

Vydac TP 系列 Avantor Hichrom 產品規格表

Phase	Particle Size (μm)	Pore Size ($\text{\AA}$)	Surface Area (m^2/g)	Carbon Loading	Endcapping	USP
218TP C18	3 5 10	300	60-110	8%	YES	L1
238TP C18	5	300	60-110	4%	YES	L1
208TP C8	5	300	60-110	5%	YES	L7
214TP C4	5 10	300	60-110	3%	YES	L26
219TP Diphenyl	5 10	300	60-110	4%	YES	L11