

HAI Speed 10 CARTRIDGES

10mm I.D. Semi-Prep Cartridge Columns



Guide to HAI Speed 10 Cartridge Part Numbers

xF-0310-xxxx 30 x 10mm cartridge
 xF-0510-xxxx 50 x 10mm cartridge
 xF-1010-xxxx 100 x 10mm cartridge

See Page 23 for complete Part Number information

HAI Speed 10 - A high performance 10mm I.D. Cartridge System

High throughput requirements for combinatorial chemistry applications often call for 10mm I.D. columns. Higgins Analytical has developed a convenient cartridge system that offers an economical solution for demanding combinatorial and semi-preparative needs. Even though a HAI Speed 10 cartridge costs a fraction of conventional 10mm HPLC column of the same length, they have the same consistently high performance of other Higgins Analytical columns.

Each 5cm and 10cm cartridge column is tested and is shipped with a CoFA.

CLIQUEUS 10mm I.D. Semi-Prep Cartridges



CLIQUEUS C18 10µm

10mm I.D. s.s. Guard Cartridges

30x10mm 2/pkg
 CF-0310-C181

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg 100x10mm 1/pkg
 CF-0510-C181 CF-1010-C181

CLIQUEUS C18 5µm

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg 100x10mm 1/pkg
 CF-0510-C185 CF-1010-C185

HAI SIL 100 10mm I.D. Semi-Prep Cartridges



HAI SIL 300 C18 10µm

10mm I.D. s.s. Guard Cartridges

30x10mm 2/pkg
 HF-0310-W181

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg 100x10mm 1/pkg
 HF-0510-W181 HF-1010-W181

HAI SIL 300 10mm I.D. Semi-Prep Cartridges



HAI SIL 300 C18 5µm

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg 100x10mm 1/pkg
 HF-0510-W185 HF-1010-W185

PROTO 200 10mm I.D. Semi-Prep Cartridges



PROTO 200 C4 10µm

10mm I.D. s.s. Guard Cartridges

30x10mm 2/pkg
 RF-0310-D041

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg 100x10mm 1/pkg
 RF-0510-D041 RF-1010-D041

PROTO 200 C4 5µm

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg 100x10mm 1/pkg
 RF-0510-D045 RF-1010-D045

PROTO 200 C18 10µm

10mm I.D. s.s. Guard Cartridges

30x10mm 2/pkg
 RF-0310-D181

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg 100x10mm 1/pkg
 RF-0510-D181 RF-1010-D181

PROTO 200 C18 5µm

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg 100x10mm 1/pkg
 RF-0510-D185 RF-1010-D185



PROTO 300 10mm I.D. Semi-Prep Cartridges



PROTO 300 C4 10µm

10mm I.D. s.s. Guard Cartridges

30x10mm 2/pkg
RF-0310-W041

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg	100x10mm 1/pkg
RF-0510-W041	RF-1010-W041

PROTO 300 C4 5µm

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg	100x10mm 1/pkg
RF-0510-W045	RF-1010-W045

PROTO 300 C18 10µm

10mm I.D. s.s. Guard Cartridges

30x10mm 2/pkg
RF-0310-W181

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg	100x10mm 1/pkg
RF-0510-W181	RF-1010-W181

PROTO 300 C18 5µm

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg	100x10mm 1/pkg
RF-0510-W185	RF-1010-W185

PHALANX 10mm I.D. Semi-Prep Cartridges



PHALANX C18 10µm

10mm I.D. s.s. Guard Cartridges

30x10mm 2/pkg
PF-0310-C181

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg	100x10mm 1/pkg
PF-0510-C181	PF-1010-C181

PHALANX C18 5µm

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg	100x10mm 1/pkg
PF-0510-C185	PF-1010-C185

TARGA 10mm I.D. Semi-Prep Cartridges



TARGA C18 10µm

10mm I.D. s.s. Guard Cartridges

30x10mm 2/pkg
TF-0310-C181

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg	100x10mm 1/pkg
TF-0510-C181	TF-1010-C181

TARGA C18 5µm

10mm I.D. s.s. Semi-Prep Cartridges

50x10mm 1/pkg	100x10mm 1/pkg
TF-0510-C185	TF-1010-C185

See pages 66 - 67 for information about Cartridge Holders required for these Cartridges

