

## HAIPEEK CARTRIDGES

### PEEK/TITANIUM CARTRIDGE COLUMNS



#### Guide to HAIPEEK Analytical Cartridge Part Numbers

xK-0421-xxxx 40 x 2.1mm PEEK/titanium cartridge  
 xK-0521-xxxx 50 x 2.1mm PEEK/titanium cartridge

See Page 23 for complete Part Number information

#### PEEK and Titanium Cartridges

The HAIPEEK family of narrowbore HPLC cartridge columns is a unique high performance product for LC-MS and other high-speed, small volume applications. The 40 and 50mm sizes of 2.1mm I.D HAIPEEK HPLC Cartridges are individually tested. A Certificate of Analysis, which includes the original chromatogram and comprehensive testing and sorbent documentation, accompanies each column. The 20mm x 2.1mm cartridges excel as narrowbore guard columns, yet even though they are not individually tested, the majority are routinely used as very economical and efficient columns for high speed HPLC screening and LC-MS applications. See pages 66 - 67 for information about the cartridge holders required for Higgins Analytical's HAIPEEK Cartridge Columns. HAIPEEK cartridges are also compatible with PE/Brownlee MPLC cartridge holders.

### CLIQUEUS 2.1mm I.D. PEEK/Titanium Cartridges



#### CLIQUEUS C8 5µm

20X2.1mm CK-0221-C085	5/pkg	40X2.1mm CK-0421-C085	2/pkg	50X2.1mm CK-0521-C085	2/pkg
--------------------------	-------	--------------------------	-------	--------------------------	-------

#### CLIQUEUS C18 5µm

20X2.1mm CK-0221-C185	5/pkg	40X2.1mm CK-0421-C185	2/pkg	50X2.1mm CK-0521-C185	2/pkg
--------------------------	-------	--------------------------	-------	--------------------------	-------

#### CLIQUEUS Cyano 5µm

20X2.1mm CK-0221-CNP5	5/pkg	40X2.1mm CK-0421-CNP5	2/pkg	50X2.1mm CK-0521-CNP5	2/pkg
--------------------------	-------	--------------------------	-------	--------------------------	-------

#### CLIQUEUS Phenyl 5µm

20X2.1mm CK-0221-PHN5	5/pkg	40X2.1mm CK-0421-PHN5	2/pkg	50X2.1mm CK-0521-PHN5	2/pkg
--------------------------	-------	--------------------------	-------	--------------------------	-------

#### CLIQUEUS Silica 5µm

20X2.1mm CK-0221-SIL5	5/pkg	40X2.1mm CK-0421-SIL5	2/pkg	50X2.1mm CK-0521-SIL5	2/pkg
--------------------------	-------	--------------------------	-------	--------------------------	-------

### HASIL 100 2.1mm I.D. PEEK/Titanium Cartridges



#### HASIL 100 C8 5µm

20X2.1mm HK-0221-M085	5/pkg	40X2.1mm HK-0421-M085	2/pkg	50X2.1mm HK-0521-M085	2/pkg
--------------------------	-------	--------------------------	-------	--------------------------	-------

#### HASIL 100 C18 5µm

20X2.1mm HK-0221-M185	5/pkg	40X2.1mm HK-0421-M185	2/pkg	50X2.1mm HK-0521-M185	2/pkg
--------------------------	-------	--------------------------	-------	--------------------------	-------



## HAIPEEK PEEK/TITANIUM CARTRIDGE COLUMNS

See pages 66 - 67 for information about the Cartridge Holders required for these Cartridges



### HAISIL 300 2.1mm I.D. PEEK/Titanium Cartridges

#### HAISIL 300 C18 5µm

20X2.1mm	5/pkg	40X2.1mm	2/pkg	50X2.1mm	2/pkg
HK-0221-W185		HK-0421-W185	\$310	HK-0521-W185	



### HAISIL HL 2.1mm I.D. PEEK/Titanium Cartridges

#### HAISIL HL C18 5µm

20X2.1mm	5/pkg	40X2.1mm	2/pkg	50X2.1mm	2/pkg
LK-0221-C185		LK-0421-C185		LK-0521-C185	



### PHALANX 2.1mm I.D. PEEK/Titanium Cartridges

#### PHALANX C18 5µm

20X2.1mm	5/pkg	40X2.1mm	2/pkg	50X2.1mm	2/pkg
PK-0221-C185		PK-0421-C185		PK-0521-C185	

### PROTO 200 2.1mm I.D. PEEK/Titanium Cartridges

#### PROTO 200 C4 5µm

20X2.1mm	5/pkg	40X2.1mm	2/pkg	50X2.1mm	2/pkg
RK-0221-D045		RK-0421-D045		RK-0521-D045	



#### PROTO 200 C18 5µm

20X2.1mm	5/pkg	40X2.1mm	2/pkg	50X2.1mm	2/pkg
RK-0221-D185		RK-0421-D185		RK-0521-D185	

200

### PROTO 300 2.1mm I.D. PEEK/Titanium Cartridges

#### PROTO 300 C4 5µm

20X2.1mm	5/pkg	40X2.1mm	2/pkg	50X2.1mm	2/pkg
RK-0221-W045		RK-0421-W045		RK-0521-W045	

#### PROTO 300 C18 5µm

20X2.1mm	5/pkg	40X2.1mm	2/pkg	50X2.1mm	2/pkg
RK-0221-W185		RK-0421-W185		RK-0521-W185	



### TARGA 2.1mm I.D. PEEK/Titanium Cartridges

#### TARGA C8 5µm

20X2.1mm	5/pkg	40X2.1mm	2/pkg	50X2.1mm	2/pkg
TK-0221-C085		TK-0421-C085		TK-0521-C085	

#### TARGA C18 5µm

20X2.1mm	5/pkg	40X2.1mm	2/pkg	50X2.1mm	2/pkg
TK-0221-C185		TK-0421-C185		TK-0521-C185	

