

Kromasil Classic Kromasil 100Å 產品規格表

Phase	Modification	Particle Size (µm)	Pore Size (Å)	Surface Area (m ² /g)	C-Content	End-capping	USP	pH Range	Temp. Max
Kromasil Classic C18	C18	1.8 2.5 3.5 5 7 10 13 16	100	320	20%	+	L1	2.0-10.0	60 °C
Kromasil Classic C8	C8	1.8 2.5 3.5 5 7 10 13 16	100	320	12%	+	L7	2.0-9.5	60 °C
Kromasil Classic C4	C4	1.8 2.5 3.5 5 7 10 13 16	100	320	12%	+	L26	2.0-9.5	60 °C
Kromasil Classic C1	C1	5	100	320	8%	+	L13	2.0-8.0	60 °C
Kromasil Classic Phenyl	Phenylbutyl	5 10 16	100	320	4.70%	+	L11	2.0-8.0	60 °C
Kromasil Classic NH2	Aminopropyl	3.5 5 7 10 13 16	100	320	14%	-	L8	2.0-9.5	60 °C
Kromasil Classic SIL	-	3.5 5 7 10 13 16	100	320	1.70%	-	L3	2.0-8.0	60 °C