

Kromasil Classic Kromasil 300Å 產品規格表

Phase	Modification	Particle Size (µm)	Pore Size (Å)	Surface Area (m2/g)	C-Content	End-capping	USP	pH Range	Temp. Max
Kromasil Classic C18	C18	5 10 16	300	110	8.70%	+	L1	2.0-9.5	60 °C
Kromasil Classic C8	C8	5 10 16	300	110	4.70%	+	L7	2.0-9.5	60 °C
Kromasil Classic C4	C4	5 10 16	300	110	2.90%	+	L26	2.0-9.5	60 °C
Kromasil Classic diC4	di-Butyl	10 16	300	110	N/A	N/A	N/A	N/A	60 °C
Kromasil Classic SIL	-	5 10 16	300	110	-	-	L3	2.0-8.0	60 °C