

Kromasil Classic Kromasil 60Å 產品規格表

Phase	Modification	Particle Size (µm)	Pore Size (Å)	Surface Area (m2/g)	C-Content	End-capping	USP	pH Range	Temp. Max
CN	Cyanopropyl	5 10 16	60	540	12%	–	L10	2.0-8.0	60 °C
Diol	Diol	5 10	60	540	10%	–	L20	2.0-8.0	60 °C
HILIC-D	Diol	5 10	60	540	10%	–	L20	2.0-8.0	60 °C
SIL	–	5 7 10 13 16	60	540	–	–	L3	2.0-8.0	60 °C