

Kromasil Eternity XT 產品規格表

Phase	Modification	Particle Size (µm)	Pore Size (Å)	Surface Area (m ² /g)	C-Content	End-capping	USP	pH Range	Temp. Max
Kromasil EternityXT C18	C18	1.8 2.5 5 10	100	330	19%	+	L1	1.0-12.0	60 °C
Kromasil EternityXT C8	C8	1.8 2.5 5 10	100	330	13%	+	L7	1.0-12.0	60 °C
Kromasil EternityXT PhenylHexyl	PhenylHexyl	1.8 2.5 5	100	330	15%	+	L11	2.0-12.0	60 °C