

Main column				Guard/Pre-column		
Code No.	Name	Column dimensions [mm]	USP	Code No.	Name	Column dimensions [mm]
Ion exchange chromatography cation exchange resin for amino acids						
0-019-01	CK10U	6×120		0-033-21	AFR2-PC	6×50
Ion exchange chromatography cation exchange resin for sugars						
0-009-01	CK08S	8×500	L58	0-009-11	CK08SG	6×50
0-010-01 0-010-06	CK08E	8×300 7.8×300	L58	0-010-11	CK08EG	6×50
0-010-02 0-010-07	CK08EC	8×300 7.8×300	L19	0-010-12	CK08ECG	6×50
Ion exchange chromatography cation exchange resin for carboxylic acids						
0-010-05 0-010-08	CK08EH	8×300 7.8×300	L17	0-010-15	CK08EHG	6×50
Ion exchange chromatography cation exchange resin for oligosaccharides						
0-001-01	CK02A	20×250	L58	0-001-11	CK02AG	8×10
0-001-02	CK02AS	20×250		0-001-12	CK02ASG	8×10
0-003-01	CK04S	10×200	L58	0-017-11 0-003-11	CK10SG CK04SG	6×50 8×10
0-003-02	CK04SS	10×200		0-017-11 0-003-12	CK10SG CK04SSG	6×50 8×10
Ion exchange chromatography anion exchange resin for carboxylic acids and sugars						
0-111-01	CA08F	4.6×250		0-111-11	CA08FG	4×10
Ion chromatography for cations						
0-034-01 0-034-04 0-407-01	SCK01 SCK01 CHK45/C05	6×50 4.6×150 4.6×150		0-034-21	SCK-PC	6×50
Ion chromatography for anions						
0-133-02	SCA04/PEEK	4.6×150	L31	0-133-12 0-130-22	SCA04G SCA-PC	4.6×30 8×10
Bioseparation for size exclusion						
0-213-01 0-214-01 0-215-01	CQP06 CQP10 CQP30	7.5×600 7.5×600 7.5×600	L25 L38 L37, 38	0-213-11 0-214-11 0-215-11	CQP06G CQP10G CQP30G	4×50 4×50 4×50

Main column				Main column			
Code No.	Name	Column dimensions [mm]	USP	Code No.	Name	Column dimensions [mm]	USP
Analytical and preparative chromatography columns for pharmaceutical applications [CHP column series]							
0-401-05	CHP20/C04	4.6X150	L21				
0-401-03	CHP20/C04	20X150	L21				
0-403-05	CHP20/C10	4.6X150	L21				
0-403-01	CHP20/C10	4.6X250	L21				
0-403-02	CHP20/C10	10X250	L21				
0-403-03	CHP20/C10	20X150	L21				
0-403-04	CHP20/C10	20X250	L21				
0-405-01	CHP07/C04	4.6X150					
0-405-04	CHP07/C04	20X200					
0-405-05	CHP07/C10	4.6X150					
0-406-01	CHP07/C10	4.6X250					
0-406-02	CHP07/C10	10X150					
0-406-03	CHP07/C10	20X150					
0-406-04	CHP07/C10	20X250					
0-402-05	CMG20/C04	4.6X150	L39				
0-402-03	CMG20/C04	20X150	L39				
0-202-06	CMG20/C10	4.6X150	L39				
0-202-05	CMG20/C10	4.6X250	L39				
0-202-02	CMG20/C10	10X250	L39				
0-202-03	CMG20/C10	20X150	L39				
0-202-04	CMG20/C10	20X250	L39				
0-404-01	CHK40/C04	4.6X150					
0-407-01	CHK45/C05	4.6X150					
0-504-01	CHPOD/C04	4.6X150					
0-504-04	CHPOD/C04	20X200					
XtalSpeed™ series							
0-047-01	DA01	4.6X50					
0-047-04	DA01	4.6X100					
0-047-02	DA01	7.5X100					
0-047-03	DA01	11.5X100					
0-047-11	SP01	4.6X50					
0-047-12	SP01	4.6X100					
0-047-13	SP01	7.5X100					
0-047-14	SP01	11.5X100					