

Ion Exchange MCI 三菱
Anion exchange resin: MCI GEL™ CA and CDR series 產品規格表

Product Name	Base Material	Functional group	Counter ion	Mean particle size (μm)	Cross linkage (%)	Ion exchange capacity meq / g	Typical Application
CA06S	St / DVB	QA	Cl ⁻	11	6	>1.2	Sugars Carboxylic acids
CA08S	St / DVB	QA	Cl ⁻	11	8	>1.2	Sugars Carboxylic acids
CA08F	St / DVB	QA	Cl ⁻	7	8	>1.2	Sugars Carboxylic acids
CA08Y	St / DVB	QA	Cl ⁻	25	8	>1.2	Sugars Carboxylic acids
CA08P	St / DVB	QA	Cl ⁻	120	8	>1.3	Sugars Carboxylic acids
CA10S	St / DVB	QA	Cl ⁻	11	10	>1.2	Sugars Carboxylic acids
CDR10	St / DVB	QA	Cl ⁻	7	-	>0.3	Nucleic acids Sugars