

ACE® HILIC Analytical Columns 產品規格表

Bonded Phase	USP Listing	Particle Size (µm)	Pore Size (Å)	pH Range	Temp. Max °C (low pH)	Temp. Max °C (high pH)	Carbon Load (%)	Surface Area (m ² /g)	Endcapped
HILIC-A	L3	1.7 3 5	100	2-7	40	40	N/A	300	No
HILIC-B	L95	1.7 3 5	100	2-7	40	40	4	300	No
HILIC-N	L95	1.7 3 5	100	2-7	40	40	7	300	No