

ACE® UltraCore BIO Solid-Core Columns 產品規格表

Bonded Phase	USP Listing	Particle Size (µm)	Pore Size (Å)	pH Range	Temp. Max °C (low pH)	Temp. Max °C (high pH)	Carbon Load (%)	Surface Area (m ² /g)	Endcapped
C18	L1	2.5 3.5	500 300	1-8	90	40	1.4 1.0	23 16	YES
C4	L26	2.5 3.5	500 300	2-9	90	40	0.6 0.4	23 16	YES
Phenyl2	L11	2.5 3.5	500 300	2-9	90	40	1.0 0.7	23 16	YES