

Partisil 系列 Avantor Hichrom 產品規格表

Phase	USP	Particle Size (µm)	Pore Size (Å)	pH Range	Temp. Max °C (low pH)	Temp. Max °C (high pH)	Carbon Load (%)	Surface Area (m2/g)	Endcapped
Silica	L27	5 10	85	2-7	40	40	N/A	350	N/A
C8	L10	5 10	85	2-8	60	40	8.5	350	YES
ODS	L1	5 10	85	2-8	60	40	5	350	No
PAC	L18	5 10	85	2-8	40	40	16	350	YES
SAX	L14	5 10	85	1.5-7.5	60	40	10	350	NO
SCX	L9	5 10	85	1.5-7.5	60	40	10	350	NO
ODS2	L1	5 10	85	2-8	60	40	16	350	NO
ODS3	L1	5 10	85	2-8	60	40	10.5	350	YES