

Partisphere 系列 Avantor Hichrom 產品規格表

Phase	USP	Particle Size (µm)	Pore Size (Å)	pH Range	Temp. Max °C (low pH)	Temp. Max °C (high pH)	Carbon Load (%)	Surface Area (m2/g)	Endcapped
Silica	L27	5	120	2-7	40	40	N/A	160	No
C8	L10	5	120	2-8	60	40	6	160	YES
C18	L1	5	120	2-8	60	40	10	160	YES
PAC	L18	5	120	2-8	60	40	N/A	160	YES
SAX	L14	5	120	1.5-7.5	60	40	N/A	160	NO
SCX	L9	5	120	1.5-7	60	40	N/A	160	NO
TAC-1 (PFP)	L43	5	120	2-8	60	40	N/A	160	YES
WAX (MAX-1)	N/A	5	120	1.5-7.5	60	40	N/A	160	NO