

Ultrasphere 系列 Avantor Hichrom 產品規格表

Phase	USP	Particle Size (µm)	Pore Size (Å)	pH Range	Temp. Max °C (low pH)	Temp. Max °C (high pH)	Carbon Load (%)	Surface Area (m2/g)	Endcapped
ODS	L1	3 5	80	2-8	60	40	12	190	YES
C18	L7	3 5	80	2-8	60	40	-	190	YES
CN	L10	3 5	80	2-8	60	40	-	190	NO
Silica	L3	3 5	80	2-7	40	40	N/A	190	N/A
C8 IP	L1	5	80	2-8	60	40	-	190	YES