

# FTIR 固體樣品測定



## KBr Spectrograde™ Powder 光譜級 KBr 粉末

Only ICL offers SpectroGrade™ crystal powders, which are ground from the purest crystal boules. SpectroGrade™ powder is made from ICL's own crystals and is free from spurious absorbances. Because ICL is the only spectroscopy accessories manufacturer that grows its own crystals, ICL can offer ultra-pure powders that will satisfy the most demanding spectroscopist. Try our new 99.999% pure FTIR and XL grades available exclusively from ICL.

Part No.	Description	Part No.	Description
0011-8966	Premium Grade FTIR KBr Powder, 100gms	0011-1391	SpectroGrade™ CsI Powder, 50 grams
0011-4540	FTIR SpectroGrade™ KBr powder, 100 grams	0011-2476	SpectroGrade™ TiBr (KRS-5) powder, 50 grams
0011-4541	XL SpectroGrade™ KBr powder 99.999% pure, 100 grams	0011-4129	SpectroGrade™ KCl powder, 100 grams
0011-4998	FTIR SpectroGrade™ KBr powder, 25 grams	0011-5078	XL SpectroGrade™ KBr powder 99.999% pure, 25 grams
		0011-4875	SpectroGrade™ NaCl Powder, 75 grams



## Random Cuttings-Chunks of KBr and KCl

### KBr 及 KCl 結晶顆粒

Keep your powder dry! And clean. How? Grind your own KBr powder from random cuttings of only the purest of ICL's crystal boules. Only the cleanest optical materials, free from spurious absorption bands, are used.

Part No.	Description	Part No.	Description
0002-2498	KBr Random Cuttings—1 pound	0002-4980	KCl Random Cuttings—50 grams
0002-4539	KBr Random Cuttings—100 grams	0002-5583	KCl Random Cuttings—25 grams
0002-5582	KBr Random Cuttings—50 grams	0002-4503	CaF2 Random Cuttings—1/2 pound
0002-2499	KCl Random Cuttings—1 pound	0002-4523	CaF2 Random Cuttings—1 pound
0002-4538	KCl random cuttings—100 grams	0002-4534	ZnS (Cleartran) Random Cuttings—1/2 pound
		0002-4789	NaCl Random Cuttings—100 grams



## Agate Mortars and Pestles 瑪瑙研鉢及杵

外徑 : 25mm, 40mm, 50mm, 60mm, 65mm,  
70mm, 75mm, 80mm, 90mm, 100mm,  
125mm, 150mm, 175mm, 200mm