

# D2燈明細表 (For UV/Vis檢測器&光譜儀) ---1---

廠 牌	機 型	D2 燈型式
ABI	120, 980, 1000S, Series 200, SF773/783/785A	AB-400
ACS	LC12/16/750	AC-400
Agilent	1050 VW (A)	HP-400
Agilent	1050 VW (C), 1050 DAD, 1100 Series, 8453	HP-404
Agilent	G1314A/B, G1315A/B, G1365A/B	HP-404
Beckman	DU 600, 6000, 620, 630, 640, 7000, 7500 , 800	BE-403
Cecil Instruments	CE 100, 1000, 1070, 2000, 5000, 6000 Series, Series 2	CE-400
Dionex	PDA 100, 3000	DI-401
Dionex/Gynkotek	340S DAD, UVD160/170S/170U/320/340S/340U/33	GY-400
GBC	LC 1200	GC-400
GBC	LC 1205/1206/1210	GC-402
Gilson	Holochrome	GI-400
Gilson	115, 116, 117, 118, 119, 151, 155, 160/170 Diode Array	GI-401
Jasco	870、875	JV-400
Jasco	915DAD、970A~、975A	JV-400
Jasco	970B~、975B~、1570、1575	JV-401
Jasco	2070 / 2075	JV-402
Knauer	2151, 8700, 9700	KA-400
Knauer	Welchrom K-2000/2001/2500/2501/2600	KA-401
Kontron	430(A)	KN-400
Kontron	332, 432, Uvikon 430/432	KN-401
Kontron	710, 810, Uvikon 600/700/800/942/943	KN-403
LDC	3000, 3200, 4000, 4200, Spectromonitor D/I/II/III	LD-400
Linear	200, 504, 206, UV100/200, Focus 2000	LI-401
Merck/Hitachi	U1000/2000/3210/4010	MH-404
Merck/Hitachi	L2500/4000/4250/7000/7400/7420/7450/7455 DAD	MH-404
Perkin Elmer	Int. 2000/4000, Lambda 1, LC135/55/75/85/90	PE-401
Perkin Elmer	LC90	PE-402

# D2燈明細表 (For UV/Vis檢測器&光譜儀) ---2---

廠 牌	機 型	D2 燈型式
Perkin Elmer	Series 200 DAD, Phase I / II	PE-403
Perkin Elmer	Lambda 2/10/25/35/45/650/850/950, LC-480 PDA	PE-404
Perkin Elmer	Lambda Bio 19/20/40/800/900	PE-404
Perkin Elmer	M55/550/552A/559A/575	PE-404
Perkin Elmer	AAAnalyst 600AA/700AA/800AA, AAS 1100/2100/4100	PE-404
Shimadzu	SPD 6A/6AV/7A, LC4A/6A	SH-400
Shimadzu	SPD10A/10AV/10AVP	SH-401
Shimadzu	AA-6200/670, CS-9000, UV-1200/1600/1700/1800/2100 series	SH-405
Spectra Physics/ Thermo Seps	UV 6000, Accela PDA, Surveyor PDA	SP-403
Spectra Physics/ Thermo Seps	SP8200/8400/8430/8440/8450/8490	SP-404
Thermo Separation	UV100/1000/150/200/2000/3000	LI-401
Varian	UV100 / 200 / 50 / 9050, Prostar 310, Vista 5000	VAR-400
Varian	Spectra A10/20/30/40, Spectra AA 880/220 FS/300	VAR-401
Varian	AA300/475/875/1275/1475/3000	VAR-401
Varian	Prostar 340/345	VAR-402
Varian	Prostar 325/335	VAR-404
Waters	480/481	WA-400
Waters	484	WA-401
Waters	486	WA-402
Waters	450	WA-403
Waters	996/2996 PDA	WA-404
Waters	2487/2488	WA-405

**廠牌型號繁多，無法一一詳列，不在表內之**

**各廠牌、型號，請來電洽本公司查詢**