



原子吸收光譜儀(AA)燈管



- * 2" 直徑燈管：適用 Perkin Elmer 各機型。
- * 1.5" 直徑燈管：適用 Varian、Hitachi、Shimadzu、GBC、Unicam...等。
- * 多元素燈管：貳元素／參元素／肆元素／伍元素／陸元素。

元素	符號	Primary	Secondary	Max	Operating	NON-DATA CODED		DATA CODED	
						1.5"	2"	1.5"	2"
Aluminum	鋁	Al	396.2nm	308.2nm	12mA	7-12mA		<input checked="" type="checkbox"/>	
Antimony	銻	Sb	217.6	206.8	15	7-10		<input checked="" type="checkbox"/>	
Arsenic	砷	As	193.7	189.0	12	7-10		<input checked="" type="checkbox"/>	
Barium	鋇	Ba	553.5	455.4	12	7-10		<input checked="" type="checkbox"/>	
Beryllium	鈹	Be	234.9	---	15	7-15		<input checked="" type="checkbox"/>	
Bismuth	鉍	Bi	223.1	222.8	10	7-10		<input checked="" type="checkbox"/>	
Boron	硼	B	249.8	208.9	15	7-15		<input checked="" type="checkbox"/>	
Cadmium	鎘	Cd	228.8	326.1	10	3-7.5		<input checked="" type="checkbox"/>	
Calcium	鈣	Ca	422.7	239.9	10	5-7.5		<input checked="" type="checkbox"/>	
Cerium	鈰	Ce	520.0	569.7	15	6-10		<input checked="" type="checkbox"/>	
Cesium	銫	Cs	852.1	455.6	10	7-10		<input checked="" type="checkbox"/>	
Chromium	鉻	Cr	357.9	425.4	15	5-10		<input checked="" type="checkbox"/>	
Cobalt	鈷	Co	240.7	304.4	15	7-15		<input checked="" type="checkbox"/>	
Copper	銅	Cu	324.8	217.9	15	5-10		<input checked="" type="checkbox"/>	
Dysprosium	鐳	Dy	419.5	404.6	15	5-10		<input checked="" type="checkbox"/>	
Erbium	鉺	Er	400.8	389.3	15	5-10		<input checked="" type="checkbox"/>	
Europium	鐳	Eu	459.4	462.7	15	7-15		<input checked="" type="checkbox"/>	
Gadolinium	釷	Gd	368.4	405.8	15	7-10		<input checked="" type="checkbox"/>	
Gallium	鎳	Ga	287.4	403.3	15	5-7.5		<input checked="" type="checkbox"/>	
Germanium	鍺	Ge	465.1	271.0	15	5-10		<input checked="" type="checkbox"/>	
Gold	金	Au	242.8	367.6	12	5-10		<input checked="" type="checkbox"/>	
Hafnium	鈦	Hf	307.8	268.2	15	5-10		<input checked="" type="checkbox"/>	
Holmium	釹	Ho	410.4	425.4	15	5-10		<input checked="" type="checkbox"/>	
Hydrogen	氫	H	Continuum		40	5-40	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>

原子吸收光譜儀(AA)燈管

---2---


- * 2" 直徑燈管：適用 Perkin Elmer 各機型。
- * 1.5" 直徑燈管：適用 Varian、Hitachi、Shimadzu、GBC、Unicam...等。
- * 多元素燈管：貳元素／參元素／肆元素／伍元素／陸元素。

元 素	符 號	Primary	Secondary	Max	Operating	NON-DATA CODED		DATA CODED	
						1.5"	2"	1.5"	2"
Indium	銦	In	303.9	325.6	15	4-7.5		<input checked="" type="checkbox"/>	
Iridium	銱	Ir	208.9	264.0	15	7-10		<input checked="" type="checkbox"/>	
Iron	鐵	Fe	248.3	248.3	15	7-15		<input checked="" type="checkbox"/>	
Lanthanum	鐳	La	550.1	403.7	15	6-10		<input checked="" type="checkbox"/>	
Lead	鉛	Pb	217.0	261.4	8	3-7.5		<input checked="" type="checkbox"/>	
Lithium	鋰	Li	670.8	323.3	10	5-7.5		<input checked="" type="checkbox"/>	
Lutetium	鐳	Lu	335.9	356.7	15	7-10		<input checked="" type="checkbox"/>	
Magnesium	鎂	Mg	285.2	202.5	10	3-7.5		<input checked="" type="checkbox"/>	
Manganese	錳	Mn	279.5	279.8	15	5-15		<input checked="" type="checkbox"/>	
Mercury	汞	Hg	253.7	---	7	3-7		<input checked="" type="checkbox"/>	
Molybdenum	鉬	Mo	313.3	320.9	15	7-10		<input checked="" type="checkbox"/>	
Neodymium	釹	Nd	492.5	463.4	15	5-10		<input checked="" type="checkbox"/>	
Nickel	鎳	Ni	232.0	231.1	15	7-15		<input checked="" type="checkbox"/>	
Niobium	鈮	Nb	334.9	405.9	15	7-15		<input checked="" type="checkbox"/>	
Osmium	銱	Os	290.9	305.9	15	7-10		<input checked="" type="checkbox"/>	
Palladium	鈀	Pd	247.6	244.8	15	3-7.5		<input checked="" type="checkbox"/>	
Phosphorous	磷	P	213.6	214.9	15	8-15		<input checked="" type="checkbox"/>	
Platinum	白金	Pt	265.9	264.7	15	7-15		<input checked="" type="checkbox"/>	
Potassium	鉀	K	766.5	404.4	10	7.5-10		<input checked="" type="checkbox"/>	
Praesodymium	鐳	Pr	495.1	513.3	15	7-10		<input checked="" type="checkbox"/>	
Rhenium	銲	Re	346.0	346.5	15	6-10		<input checked="" type="checkbox"/>	
Rhodium	銲	Rh	343.5	328.1	15	5-10		<input checked="" type="checkbox"/>	
Rubidium	銣	Rb	780.0	794.8	12	6-10		<input checked="" type="checkbox"/>	
Ruthenium	鈳	Ru	349.9	392.6	15	5-10		<input checked="" type="checkbox"/>	
Samarium	釷	Sm	429.7	476.0	15	7-15		<input checked="" type="checkbox"/>	
Scandium	釷	Sc	391.2	390.8	15	5-10		<input checked="" type="checkbox"/>	
Selenium	硒	Se	196.0	204.0	7	5-7		<input checked="" type="checkbox"/>	
Silicon	矽	Si	251.6	250.7	15	7-15		<input checked="" type="checkbox"/>	
Silver	銀	Ag	328.1	338.3	10	5-10		<input checked="" type="checkbox"/>	
Sodium	鈉	Na	589.0	330.2	10	5-7.5		<input checked="" type="checkbox"/>	
Strontium	銣	Sr	460.7	407.8	12	7-10		<input checked="" type="checkbox"/>	
Tantalum	鉭	Ta	271.5	275.8	15	5-10		<input checked="" type="checkbox"/>	
Tellurium	碲	Te	214.3	225.9	10	7-10		<input checked="" type="checkbox"/>	

- * 2" 直徑燈管：適用 Perkin Elmer 各機型。
- * 1.5" 直徑燈管：適用 Varian、Hitachi、Shimadzu、GBC、Unicam...等。
- * 多元素燈管：兩元素／參元素／肆元素／伍元素／陸元素。

元素	符號	Primary	Secondary	Max	Operating	NON-DATA CODED		DATA CODED	
						1.5"	2"	1.5"	2"
Terbium	鈿	Tb	432.7	431.9	15	5-10		<input checked="" type="checkbox"/>	
Thallium	鉍	Tl	276.7	258.0	15	6-10		<input checked="" type="checkbox"/>	
Thorium	鈷	Th	371.9	---	20	6-10		<input checked="" type="checkbox"/>	
Thullium	錕	Tm	371.8	436.0	15	6-10		<input checked="" type="checkbox"/>	
Tin	錫	Sn	233.5	224.6	15	5-7.5		<input checked="" type="checkbox"/>	
Titanium	鈦	Ti	364.3	365.4	15	7-10		<input checked="" type="checkbox"/>	
Tungsten	鎢	W	255.1	294.7	15	10-15		<input checked="" type="checkbox"/>	
Uranium	鈾	U	358.5	356.6	15	6-10		<input checked="" type="checkbox"/>	
Vanadium	釩	V	318.5	306.6	15	7-10		<input checked="" type="checkbox"/>	
Ytterbium	鐿	Yb	398.8	346.4	12	5-7.5		<input checked="" type="checkbox"/>	
Yttrium	釷	Y	410.2	414.2	10	5-7.5		<input checked="" type="checkbox"/>	
Zinc	鋅	Zn	213.9	307.6	15	4-10		<input checked="" type="checkbox"/>	
Zirconium	鋯	Zr	360.1	468.7	15	7-15		<input checked="" type="checkbox"/>	

Multi-Element Lamps (多元素燈管)



多 元 素	
Ca/Mg	Cu/Cr/Ni/Ag
Ca/Zn	Cu/Cr/Fe/Ag
K/Na	Cu/Fe/Mn/Zn
Ag/Au	Cu/Co/Cr/Mn/Ni
Al/Ca/Mg	Cu/Co/Fe/Mn/Mo
Ca/Mg/Zn	Cu/Co/Cr/Ni/Ag
Cu/Fe/Ni	Cu/Co/Cr/Fe/Mn/Ni

Adapters and Power Supply (連接器&電源供應器)

